



**DATA SHEET** 

**Model: ESLE-8000**

**Binocular head loupes with lamp**



- Diameter of the light spot at a distance of 420mm: 20-90mm
- Illumination in 250mm working distance: >30000 lx
- Fluctuation distance: 12.5mm
- Vertical direction of Illumination: +/-25 degree adjustable
- Net Weight 70g (excluding headband)
- Illumination Source: extremely bright light, 5W LED. Life 20000 hours continuous
- Color temperature of 6000 Kelvin
- LED lamp without the heat of traditional Halogen headlight
- Smart design, convenient for mobile diagnosis.
- Magnification: 3.5x
- Working distance: 420 mm
- Field of view: 42 mm
- Battery: 7.2V Lithium battery (rechargeable)
- The headlight could be used when it is charging.
- Low-battery indicator light.
- Battery Run-time: >2.5 hours
- Voltage: 110VAC/ 220 VAC
- Maximum Power input: 15W

**CE**



**ESSE3** srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

