




# VISUAL ACUITY CHART



Cod. **ESWH-1401**

- Reading Distance: 5 m
- Power supply: 220V 50Hz
- Dimensions: 800 x 450 x 18 mm



**MULTIFUNCTION VISION CHART**  
TO BE TESTED AT 5 METERS (16.4 FEET)

5-LogMAR DECIMAL	20/FOOT 6/METER	Visual Acuity Test Patterns
4.0 0.1	20/200 6/60	Large 'E' character; Red and green color discrimination test (8, 9, 6, 3, 2, 0, 8, 9, 6, 0, 2, 3)
4.1 0.12	20/160 6/48	'M' and '3' characters; Clock face test
4.2 0.15	20/125 6/38	'E' and 'W' characters; Grid of dots test
4.3 0.2	20/100 6/30	'E', 'W', and '3' characters; Diamond and cross test
4.4 0.25	20/80 6/24	'3', 'M', and 'W' characters; Red dot test
4.5 0.3	20/63 6/19	'W', '3', 'M', and 'E' characters
4.6 0.4	20/50 6/15	'E', 'M', '3', and 'W' characters
4.7 0.5	20/40 6/12	'E', 'W', '3', and 'M' characters
4.8 0.6	20/32 6/9.5	'3', 'W', 'M', 'E', 'W', '3' characters
4.9 0.8	20/25 6/7.5	'E', 'M', 'W', 'M', 'E', 'W', '3' characters
5.0 1.0	20/20 6/6	'3', 'M', 'W', 'E', 'M', '3', 'W', 'E' characters
5.1 1.2	20/16 6/4.8	'E', 'W', '3', 'W', 'M', 'E', '3', 'M' characters
5.2 1.5	20/12.5 6/3.8	'W', '3', 'M', 'E', 'W', 'E', 'M', '3' characters




**ESSE3** srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)






# ECHELLE D'ACUITE VISUELLE



Cod. **ESWH-1401**

- Distance de lecture : 5 m
- Alimentation : 220V 50Hz
- Dimensions: 800 x 450 x 18 mm



**ESSE3** srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

