



DATA SHEET OPTICAL RETINOSCOPY RACK LENS

Cod. **RS-6R30P** Optical Retinoscopy Rack Lens

Sphere: 15pcs each of concave and convex.
0.5 1.00 1.50 2.00 2.50 3.00 3.50
4.00 4.50 5.00 6.00 7.00 8.00 9.00 10.00



Cod. **RS-4R28W** Optical Retinoscopy Rack Lens

0.50D 1.00D 1.50D 2.00D 2.50D 3.00D 3.50D
4.00D 5.00D 6.00D 7.00D 8.00D 9.00D 12.00D



Cod. **RS-4R28AL** Optical Retinoscopy Rack Lens

Aluminium case , and aluminium frame.
0.50D, 1.00D 2.00D 2.50D 3.00D 3.50D
4.00D 5.00D 6.00D 7.00D 8.00D 9.00D 12.00D



Cod. **RS-2R32PL** Retinoscopy Rack set

Soft case
0.50D 1.00D 1.50D 2.00D 2.50D 3.00D 3.50D 4.00D
....
5.00D 6.00D 7.00D 8.00D 9.00D 10.00D 12.00D 15.00D



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

