



ESCS-2000A

A Scan



OPHTHALMOLOGY / ENT



Technical specification

Ophthalmic A scan Ophthalmic biometer scan

Probe: 10MHz

Storage: 8 groups

Accuracy: 0.05mm

IoL Calculation: SRK-II, SRK/T Hoffer-Q, Holladay, SCDK, Haigis etc

Measure: ACD, Lens, Vitreous, Axial Length, and its Average

Modes: Manual, Automatic

Printer: Inner therma Printer

Power: AC220V± 10%50Hz

Display Monitor: High speed LCD

Report Data: Patient's name, ID,

Right/Left, Doctor's name, Date, IOL Table



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3.dreamgest.com

