



**DATA SHEET** 

**Model: HSL-900**

**Portable slit lamp**



#### Features and advantages

Images real-time display, 10 Times magnification

High definition imaging image resolution , 16MP

Micro SD memory card up to 16GB

USB and WiFiconnection

Low weight 600g , movable and portable

Rechargeable battery provides up to more than 4 hours of continuous operation.

3.97 inch TFT-lcd color and touch screen

Auto focus ,White LED illumination

With built-in Android APP inside . do not need other mobile or laptop

With software , can install on laptop, slit lamp connect to mobile or laptop wifi

#### Performance and Parameters

Magnification: 10X

Refractive compensation: -7D to +7D

Working distance: 80mm

Slit width: 0-10mm continuously adjustable

Aperture diameter: 1mm,5mm, 10mm, 14mm

Filter: Heat-absorption,Gray, Red-free, Cobalt blue

Light source: White LED

Image resolution: 16MP

Focus Mode: Auto

Screen: 3.97" color and touch screen

Storage: 16GB MicroSD card

Power supply: 3.7VLithium Battery 3400mAh

Interface: Mini USB/Wifi/BT/4G

N. Weight: 600G

G. Weight: 3kg

Packing size: 400mm\*200mm\*230mm



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

