



# ESLM-1

Slit lamp



OPHTHALMOLOGY / ENT



# ESLM-1

## Slit lamp

### Specification

**Type:** Converging Stereoscope  
**Total Magnification:** 10X: 16X (10X eyepiece) - 16X: 25.6X (16X eyepiece)  
**Slit Width:** 0mm-10mm continuously adjustable  
**Slit Height:** 1mm-10mm continuously adjustable  
**Aperture Diameter:**  $\Phi$  10,  $\Phi$  8,  $\Phi$  5,  $\Phi$  3,  $\Phi$  2,  $\Phi$  1,  $\Phi$  0.2mm  
**Slit Angle:** 0° -180°  
**Illumination Tilting:** 0° , 5° , 10° , 15° , 20°  
**Illumination Bulb:** Halogen bulb 12V/50W  
**Diopter Adjustment:** +5D--5D  
**Filter:** Heat absorption, grey, redfree(green), cobalt blue  
**Input Voltage:** 110V/220V 60/50Hz  
**Fixation:** Red LED  
**Power:** 60VA  
**Optional:** LED lamp, Iphone or Sumsung mobile adapter, Electric table.

### Optional



Mobile adapter



Electric table



ESSE3 srl, Via Garibaldi 30  
 14022 Castelnuovo D.B. (AT)  
 Tel +39 011 99 27 706  
 Fax +39 011 99 27 506  
 e-mail esse3@chierinet.it  
 web: www.esse3-medical.com

