



ESLM-2

Slit lamp



OPHTHALMOLOGY / ENT



ESLM-2

Slit lamp

Specification

Type: Parallel Galilean
Eyepiece: 12.5X
Total Magnification: 6X, 10X, 16X, 25X, 40X
Diameter of Visual Field: Φ 37, Φ 23, Φ 14, Φ 8.7, Φ 5.7mm
Diopter Adjustment: +5D--5D
Slit Width: 0mm-14mm continuously adjustable
Slit Height: 1mm-14mm continuously adjustable
Aperture Diameter: Φ 14, Φ 10, Φ 5, Φ 3, Φ 2, Φ 1, Φ 0.2mm
Slit Angle: 0°-180°
Illumination Tilting: 0°, 5°, 10°, 15°, 20°
Filter: Heat absorption, grey, redfree(green), cobalt blue
Lamp: German OSRAM Halogen Tungsten Lamp
Longitudinal movement: 90mm
Lateral movement: 100mm
Vertical movement: 30mm
Fine base movement: 15mm

Input Voltage: 110V/220V 60/50Hz
Input Power: 60VA
Illumination Bulb: 12V/50W halogen bulb ,LED lamp optional
Fixation: Red LED

Optional



Mobile adapter



Electric table



ESSE3 srl, Via Garibaldi 30
 14022 Castelnuovo D.B. (AT)
 Tel +39 011 99 27 706
 Fax +39 011 99 27 506
 e-mail esse3@chierinet.it
 web: www.esse3-medical.com

