



# WATER BATH

Model	Inner dimensions	Overall dimensions	weight	capacity	consumption
<b>BS/1DG</b>	30 x 16 x 16 cm	47 x 22 x 24H cm	5 kg	<b>8 lt</b>	500 W
<b>BS/2DG</b>	30 x 30 x 16 cm	53 x 38 x 26H cm	9 kg	<b>15 lt</b>	900 W
<b>BS/3DG</b>	30 x 45 x 16 cm	63 x 37 x 29H cm	10 kg	<b>28 lt</b>	1500 W
<b>BS/6DG</b>	30 x 90 x 16 cm	112 x 38 x 26H cm	17 kg	<b>45 lt</b>	2300 W

- Entirely frame in stainless steel
- Tank capacity 8 - 15 - 28 - 36 and 45 litres , heavy thermal insulation
- Electronic PID adjustable temperature up to 99 °C , precision  $\pm 0,1$  °C
- Digital display temperature, pre-setting temperature value
- Drain system with stainless steel lid ( *optional for BS/1DG and BS/2DG* )
- Power supply 220 V 50 Hz

#### Optional :

- Stainless steel cover
- Test tube holder dimensions 28 x 14 cm , 50 / 60 / 78 places diameter 13 mm or 17 mm
- Safety thermostat
- Water mixer for temperature homogenisation
- External glass thermometer