



# 89BH5@5-F 7CA DF 9GGCF G



DENTISTRY

# DENTAL AIR COMPRESSORS



## Cod. ES-30

Air flow	70L/min
Adjustable Pressure	2 to 8 Bar
Noise	53 dB
Tank capacity	30L
Power consumption	550W - 2,4V
Power supply	220V AC,50 HZ Optional: 220V AC,60HZ or 110V AC,60HZ
Dimensions	45x45x63 cm
Weight	30kg



## Cod. ES-35

Air flow	100L/min
Adjustable Pressure	2 to 8 Bar
Noise	53 dB
Tank capacity	35L
Power consumption	840W – 3,7A
Power supply	220V AC,50 HZ Optional: 220V AC,60HZ or 110V AC,60HZ
Dimensions	52x52x75 cm
Weight	35kg



## Cod. ES-45

Air flow	140L/min
Adjustable Pressure	2 to 8 Bar
Noise	53 dB
Tank capacity	45L
Power consumption	1100W – 4,4A
Power supply	220V AC,50 HZ Optional: 220V AC,60HZ or 110V AC,60HZ
Dimensions	45x90x70 cm
Weight	55kg