

# INTRAORAL CAMERA

## item: 582/00 - MD8103O INTRAORAL CAMERA



### TECHNICAL FEATURES:

Image sensor:	1/4" sony CCD
Image resolution:	2.0 mega pixels
Output:	Video/S-Video/VGA
Lamp:	6 LED (5600 K)
Accumulation point:	3 mm - 50 mm
Visual angle:	105°
Power supply:	1600mA
Frequency:	5.8 GHz
Wireless distance:	10m
Signal format:	NTS/PAL (opzionale)

### Configuration System



## item: 577/00 - MD8503O - INTRAORAL CAMERA WIRE 2.0

### TECHNICAL FEATURES:

Image sensor:	1/4" sony CCD
Image resolution:	2.0 mega pixels
Output:	USB/VGA/Uscita Video
Lamp:	6 LED (5600 K)
Accumulation point:	3 mm - 50 mm
Visual angle:	105°
Memory images:	28 in high resolution or 112 in low resolution
Cable length:	2 m
Signal format:	NTS/PAL (opzionale)
Power supply:	DC 5V 1200mA

### Configuration System



## ARMS & SUPPORTS



### MONITOR TFT LCD ARM

**item 87010000**

Monitor arm (L.385 mm), SVGA, video and power supply cable.

**item 87010000/NC**

Monitor arm (L.385 mm).

**item 87030000**

Short Monitor arm (L.130 mm), SVGA, video and power supply cable.

**item 87030000/NC**

Short Monitor arm (L.130 mm).

### CABLE

**item 87010023**

Complete cable for monitor arm 5 mt.

**item 87010028**

Complete cable for monitor arm 10 mt.

### CLAMP

**item 87010013**

Clamp for monitor arm Ø 45 mm.

**item 87010014**

Clamp for monitor arm Ø 50 mm.

**item 87010022**

Clamp for monitor arm Ø 55 mm.

**item 87010036**

Clamp for monitor arm Ø 60 mm.

**item 87010043**

Rubber clamp shim.

### ECONOMIC MONITOR ARM

**item 382/00**

Clamp for monitor arm Ø 45 mm.

### MONITOR NEOVO LCD 17"

**item 870500016**

Monitor Neovo X-17 no medical certified.

**item 870500017**

Monitor Neovo DR-17 medical certified.

## DENTISTRY MICROMOTORS



CE0476

**item: 551/00 - 4000 COLIBRÌ INTRAMATIC**

Dentistry micromotor with intramatic connection.

**item: 553/00 - 4000 COLIBRÌ INTRAMATIC R**

#### TECHNICAL FEATURES:

Power:	100 Watt
Speed:	1.000-30.000 rpm
Feeding:	230 V - Hz 50/60
Weight:	3,5 Kg