

# GYNEACOLOGICAL / DELIVERY BEDS



## ES-045 GYNECOLOGICAL EXAMINATION BED

This bed is designed for gynecological diagnosis and parturition.

It can be folded and easily dismantled.

Specifications:1870×520×740mm

Back Section Raised from the Horizontal:≥55°

Back Section Lowered from the Horizontal: ≥20°

Seat Section Raised from the Horizontal: ≥20°

Leg Section Lowered from the Horizontal: ≥85

ES-045



## ES-045S GYNECOLOGICAL EXAMINATION BED

Simplified version



## ES-046 ORDINARY DELIVERY BED

This bed is designed for obstetrics diagnosis , parturition , , dilatation and curettage as well as for gynecological operations.

Specifications:1900×600×800mm

Back Section Raised from the Horizontal:≥45°

Back Section Lowered from the Horizontal: ≥ 7°

Seat Section Raised from the Horizontal: ≥25°

Seat Section Lowered from the Horizontal: ≥8°

Leg Section Raised from the Horizontal: ≥25°

Leg Section Lowered from the Horizontal: ≥25°

ES-046



## ES-078 MULTI-PURPOSE DELIVERY BED

This bed is designed for parturition and gynecological operation.

It include two parts :the bed and the legs. The bed can be raised or lowered hydraulically by means of an oil pump.

The ascending and descending movements of the legs are controlled by the leading screw and operated by the hand wheel.

Specifications:1900×600×800(1000)mm

Back Section Raised from the Horizontal: ≥35°

Back Section Lowered from the Horizontal: ≥10°

Seat Section Raised from the Horizontal: ≥80mm

ES-078



## ES-300 MULTI-PURPOSE OPERATING / DELIVERING TABLE

Side-controlled

This table is suitable for surgical procedures of the head., neck , thorax and abdomen, perineum and extremities as well as for gynecological, ENT and orthopedics operation.

It can be raised or lowered hydraulically by means of an oil pump.

Specifications:2100×480×800-1040mm

Trendelenburg : ≥20°

Reversed Trendelenburg : ≥30°

Lateral Tilt: ≥18°

Back Section Raised from the Horizontal: ≥75°

Back Section Lowered from the Horizontal: ≥18°

Leg Section Lowered from the Horizontal: ≥90°

Kidney Bridge Elevation: ≥120mm

ES-300



ESSE 3 Via Garibaldi 30  
14022Castelnuovo D.B. (AT)  
tel +39 011 99 27 706  
fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web : [www.4ci.com](http://www.4ci.com)

