



ES-10SA

Single Flow



OXYGEN CONCENTRATOR

Model	ES-10SA (with Atomization)
Flow rate	0-10L/min
Purity	93 ± 3%
Outlet Pressure (Mpa)	0.04 - 0.07
Sound Level	≤ 50db
Power	AC220V±10%, 50/60±1Hz; AC110V±10%, 50/60±1Hz
Power Consumption	≤ 480W
Timer	Present Working Time, Accumulating Time, Presetting Time from 10mins to 40hours
Alarm	Power failure alarm, High & Low pressure alarm
Net Weight	27Kg
Size (mm)	365×375×600mm

OPTIONAL

*Low purity alarm: when oxygen purity is above 82%, it will give green light; when the purity is below 82% (82% not included), it will give red light.

*Maintenance alarm: alarm for maintenance after 3000 hours operation.

*Alarm for high temperature inside systems above 50°C .

*Pulse oximeter: fingertip pulse oximeter is available for monitoring patient oxygen saturation on site.

*Small LCD display: accumulating time; present working time presetting time.

*Dual flows.



ES-10DA

Double Flows



OXYGEN CONCENTRATOR

Model	OXY-10DA (Atomization)
Flow rate	0-10L/min
Purity	93 ± 3%
Outlet Pressure (Mpa)	0.04 - 0.07
Sound Level	≤ 50db
Power	AC220V±10%, 50/60±1Hz; AC110V±10%, 50/60±1Hz
Power Consumption	≤ 480W
Timer	Present Working Time, Accumulating Time, Presetting Time from 10mins to 40hours
Alarm	Power failure alarm, High & Low pressure alarm
Net Weight	27Kg
Size (mm)	365×375×600mm

OPTIONAL

*Low purity alarm: when oxygen purity is above 82%, it will give green light; when the purity is blow 82%(82% not included), it will give red light.

*Maintenance alarm: alarm for maintenance after 3000 hours operation.

*Alarm for high temperature inside systems above 50°C.

*Pulse oximeter: fingertip pulse oximeter is available for monitoring patient oxygen saturation on site.

*Small LCD display: accumulating time; present working time presetting time.



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.esse3-medical.com

