



ES-AUN SERIE

N series - Vacuum Steam Autoclave



3 New Technologies

- * Automatic and safety door lock system**
- * LCD display, easy for operation**
- * Steam compensation techniques**



CLASS N Model - sterilize solid instruments

- * LCD display. 2 programs, easy to choose
- * Automatically feeding water, excess water and air to be drained away
- * Display the failure reason on LCD screen. Easy to correct
- * Safety Devices:
 - a. Machine has safety door lock system
 - b. Over heat control: when inner temperature is over heated, power will be switch off automatically and start alarm system
 - c. Over pressure control: alarm will be announced when the machine is over pressure
 - d. Low water level detected: alarm will be announced when it is low water

Model	ES-14AUN	ES-17AUN	ES-23AUN
Chamber Dimension (mm) (Φ x D)	25 x 28	25 x 38	25 x 49
Chamber Capacity (liter)	14L	17L	23L
Overall Size (cm) (L x W x H)	48 x 48 x 42	48 x 53 x 42	48 x 69 x 42
Heater	1200W	1200W	1500W
Net Weight (KGS)	31	34.5	41
Gross Weight (KGS)	38	43	50
Packing Size (cm)	68 x 60 x 50	73 x 60 x 50	87 x 60 x 50
Minimum Table Size (cm)	48 x 35	48 x 40	48 x 50
Power Supply	230V, 110V, 126V		
Sterilize Temperature	121 °C / 134 °C		
Sterilize Time	Please refer to different programs		
Sterilizer Construction	Stainless Steel SUS#304		
Stainless steel instrument contain	2 trays		



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it

web www.esse3medical.com

