

ES-EX BENCH TOP STEAM STERILIZER



**ECONOMICAL
SEMI-AUTOMATIC
CONTROLLER SYSTEM**



ES-EX BENCH TOP STEAM STERILIZER

ECONOMICAL

FEATURES:

1. HORIZONTAL SEMI-AUTOMATIC STERILIZER.
2. WITH DRY AND RE-DRY FUNCTION, AUTO EXHAUST.
3. WITH LOW WATER INDICATOR.
4. WITH AUTO EXHAUST FUNCTION
5. WITH STERILIZATION TEMPERATURE 121C AND 134C BOTH.
6. DOOR CLOSE KNOB IS 90 DEGREE CLOSE/OPEN DESIGN: LIGHT AND CONVENIENT FOR OPERATION.

MODEL	ES-16EX	ES-25EX	ES-40EX
OVERALL	525W x 524D x 425H MM	525W x 691D x 425H MM	575W x 795D x 460H MM
CHAMBER DIAMETER	261MM	261MM	315MM
CHAMBER LENGTH	300MM	470MM	520MM
CHAMBER CAPACITY	16 LITERS	25 LITERS	40 LITERS
DESIGN TEMPERATURE	140 C	140 C	140 C
WORKING STERILIZATION TEMPERATURE/PRESSURE	121 °C (1,054 kg/cm2) 134 °C (2,07 kg/cm2)	121 °C (1,054 kg/cm2) 134 °C (2,07 kg/cm2)	121 °C (1,054 kg/cm2) 134 °C (2,07 kg/cm2)
FILLING WATER	MANUAL	MANUAL	MANUAL
DRY FUNCTION	YES	YES	YES
RE-DRY FUNCTION	YES	YES	YES
AUTO EXHAUST	YES	YES	YES
POWER CONSUMPTION	1.8KW	2.1KW	3.0KW
90 DEGREE DOOR OPEN/CLOSE KNOB	YES	YES	YES
OVER TEMPERATURE PROTECTION	YES	YES	YES
OVER PRESSURE PROTECTION	YES	YES	YES
ELECTRIC OVERLOAD PROTECTION	NO-FUSE BREAKER	NO-FUSE BREAKER	NO-FUSE BREAKER
CHAMBER WATER LEVEL SENSOR	YES	YES	YES

STANDARD ACCESSORIES:

1. SUS#304 STAINLESS STEEL STERILIZATION BOX x 1
2. SUS#304 STAINLESS STEEL HEATER COVER x 1

OPTIONAL ACCESSORIES:

1. SUS#304 STAINLESS STEEL STERILIZATION BASKET OR
2. TRAY SET (INCLUDING: RACK(SUS304) x1 + TRAY (ALUMINUM) x 3 + PULLING HOLDER(SUS304) x 1)



ESSE 3 Via Garibaldi 30
 14022Castelnuovo D.B. (AT)
 tel +39 011 99 27 706
 fax +39 011 99 27 506
 e-mail esse3@chierinet.it
 web : www.esse3-medical.com

