




PORTABLE OXYGEN FILLING COMPRESSOR



Cod. ES-SC



A highly efficient, portable oxygen filling Compressor.

It works with Oxygen Generator and it can be used for filling various kinds of oxygen bottles. Perfect for hospital, clinical and nursing centers.

TECHNICAL SPECIFICATIONS	
Working temperature	5°C to 35°C
Working humidity range	15% to 95% RH (condensation)
Storage temperature range	40°C to 70°C
Store humidity range	Between 10% and 95% RH (condensation)
Working voltage range	220 ± 10% V ~ 50 hz
Power range	< 600W
Insulation class	class B
Input flow/pressure	3 ~ 10 LPM @ 40 ~ 100 kpa
Bottle filling pressure	0 psi ~ 2000 ± 200 psi
Weight	About 27 kg
Filling medium	Oxygen, dry air, carbon dioxide and nitrogen
Filling time	Portable cylinder filling psi from 0 to 2000 ± 200 Psi: M6 (water volume is about 1.2 L)... 30 minutes




ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com





COMPRESSEUR PORTABLE DE REMPLISSAGE D'OXYGENE



Cod. ES-SC



Compresseur portable de remplissage d'oxygène à rendement élevé.
Il fonctionne avec le générateur d'oxygène et il peut être utilisé pour le remplissage de différentes bouteilles d'oxygène. Parfait pour les hôpitaux, cliniques et centres de soins infirmiers.

TECHNICAL SPECIFICATIONS	
Température de travail	5°C à 35°C
Taux d'humidité de travail	15% à 95% RH (condensation)
Plage de température de stockage	40°C à 70°C
Plage d'humidité de stockage	Entre 10% et 95% RH (condensation)
Plage de tension de travail	220 ± 10% V ~ 50 hz
Puissance	< 600W
Classe d'isolement	classe B
Débit /pression d'entrée	3 ~ 10 LPM @ 40 ~ 100 kpa
Pression de remplissage du cylindre	0 psi ~ 2000 ± 200 psi
Poids	Environ 27 kg
Remplissage moyen	Oxygène, air sec, dioxyde de carbone et azote
Temps de remplissage	Remplissage du cylinder: psi 0 à 2000 ± 200 Psi: M6 (le volume de l'eau est environ 1.2 L)... 30 minutes



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

