

## SUCTION UNIT-A

9GIGI 7%



### Detailed Product Description

1. The suction systems with condensate separator, air by-pass valve and electrical control box.
2. Optional power/performance ratio of the dry suction systems.
3. Responsible, energy-saving, quiet and clean vacuum supply
4. Exalted efficiency by reduction of the water consumption and the operating expenses

### Product Description:

The surgical device is more ergonomic to locate as close as possible to the dental unit and to the left side of the dental chair.

It's metal cabinet is designed for maximum ergonomym and minimum acustic and noise level. The cabinet is produced for maksimum air circulation against heating by the vacuum motor.

Cabinet surface is painted electrostatic epoxy and has smooth, bright and hygenic surface. Also the other parts are galvanized for corsion free.

The surgical aspirator (dental high vacuum system) has anti vibration system and silencer for the exhauster air and non stop suction motor's noise.

Voltage	110V/ 220V/ 230V/ 240V/380V
Motor Power	400W
Vacuum Pressure	200 mmHg
Rated Current	3 Ampers
Air flow	1000 L/min
Motor RPM	2800 RPM
Packing size	64*46*70cm
Unit size	60*42*66.5cm
G.W	40KG
N.W	35KG



ESSE 3 Via Garibaldi 30  
14022Castelnuovo D.B. (AT)  
tel +39 011 99 27 706  
fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web : [www.4ci.com](http://www.4ci.com)



## SUCTION UNIT-B



### Detailed Product Description

Advantages:

1. The suction systems with condensate separator, air by-pass valve and electrical control box.
2. Optional power/performance ratio of the dry suction systems.
3. Responsible, energy-saving, quiet and clean vacuum supply.
4. Exalted efficiency by reduction of the water consumption and the operating expenses.

### Product Description:

- 1, The system is proved to be capable of generating enough suction after strict testing.
- 2, The tubing and connecting are tightly sealed and simple in structure which make it easy to install and uninstall.
- 3, Fully Automatic control is realized. The product can start, control and turn off the drainage of the sullage automatically.
- 4, The product simultaneously collects and disposes the sullage and drains out for environment-friendly pupose.
- 5, The product is designed to turn off postponedly after the usage to suck out oral secretion, blood, tooth debris and amalgam, etc.
- 6, The product is small size, light weight, silent and easy to install.

Voltage	110V/ 220V/ 230V/ 240V/380V
Motor power	400W
Vacuum Pressure	200 mmHg
Rated Current	3 Ampers
Air flow	1000 L/min
Motor RPM	2800 RPM
Packing size	47*44*48
Unit size	43X40X44cm
G.W	29KG
N.W	24KG