

# ES-ETP



## ES-ETP

Stretcher for specialistic examination and tests.

Completely made of stainless steel.

Mobile on 4 antistatic castors provided with selector for directionality and locking system.

Radiotranslucent top equipped with guide rails that allow the insertion of X-Rays film cassette under all the top length.

Available with **manual ( standard )** or **electric control**.

U Base with lateral column.

Up/Down 660-960 mm

Trendelenburg / Reverse Trendelenburg 20°

Back section adjustment +75°

Table top size : Large 590 mm / Length 2100 mm

Mattresses are removeables and sterilizeables in autoclave

Table tops is available divided in 2 or 3 sections

Table top completely X-Ray translucent

Option

Electric control – On demand

Lateral tilting - On demand

Legs section adjustment - On demand



ESSE 3 Via Garibaldi 30  
14022Castelnuovo D.B. (AT)  
tel +39 011 99 27 706  
fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web : [www.4ci.com](http://www.4ci.com)

