



# BENCH TOP CYTO CENTRIFUGE

## ESC-CYTO CENTRIFUGE

MODEL	ESC CYTO CENTRIFUGE
MAX SPEED	3000 rpm
MAX RCF	730 xg
AMOUNT OF SINGLE TUBE	0,1 – 0,5 ml
COUNT OF ALL TRUBES	6
NOISE LEVEL	< 65 dB (A)
POWER SUPPLY	220V 50/60 Hz 3A
DIMENSIONS	300x280x270H mm
WEIGHT	15 kg



Microprocessor control, digital display , touch panel, DC brushless motor without carbon dust pollution, and many kinds of brackets for convenient use

- ◆ Electronic door interlock for safe operation.
- ◆ RCF value can be calculated automatically and showed during the run
- ◆ Self-diagnosis system can automatically detect the faults of overspeed, over temperature , imbalance and door lid interlock with displaying of the error codes on the display screen to improve the safety performances



ESSE 3 Via Garibaldi 30  
14022Castelnuovo D.B. (AT)  
tel +39 011 99 27 706  
fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web : [www.esse3-medical.com](http://www.esse3-medical.com)

