



BLOOD TYPE CARD CENTRIFUGE

ESCT-D24K



Features:

Superior performance, safety and reliability

- ◆ DC brushless motor, Microprocessor control.
- ◆ Specialized RCF button to show RCF and speed, the parameters can be revised freely during a run. It can reach the best centrifugal effect.
- ◆ Low noise and no carbon dust pollution.
- ◆ The function of electronic door interlock, over speed and self-examination for the safety of people and instrument.
- ◆ With touch panel for easier operation.
- ◆ The standard programs are build-in for blood type determination and all kinds experiment of serology, which standardize the test and research of clinical. The programs also can be set by users to satisfy the actual requirement

MODEL	ESCT-D24K
MAX SPEED	4000 rpm
Centrifugal speed of blood type card	3840 rpm
MAX RCF	1790xg
NUMBER OF MICRO GEL (SIX COLUMN)	24
SPEED CONTROL ACCURACY	± 30 rpm
TIME SETTING RATE	0 to 99min
BRUSHLESS MOTOR	AC FREQUENCY CONVERSION
NOISE LEVEL	< 60 dB (A)
POWER SUPPLY	220V 50/60 Hz
DIMENSIONS	370x360x380H mm
WEIGHT	16 Kg



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chiermet.it
web : www.esse3-medical.com

