



# ES8-\$) B

Foldable dental chair



DENTISTRY



# EG8!\$) 6

Foldable dental chair



Cod. **ESD-05B**

- High-strength aluminium alloy structure, light and stable.
- Foldable for saving space and carry it around easily.
- Comfortable and safe for up to 200kg.
- It also can be used as standard dental chair.
- Seat adjustable height: from 35 to 60 cm.
- Dimensions: 550x260x900mm.
- Weight: 7 kg



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

