


ESG-02

GYNECOLOGICAL / OBSTETRIC TABLE

FEATURES:

1. FULLY HYDRAULIC CONTROL SYSTEM.
2. THIS IS A MULTI-PURPOSE TABLE. WITH SOME OPTIONAL ACCESSORIES, IT CAN BE APPLIED AS A GYNECOLOGY EXAM, OBSTETRIC, OR PROCTOLOGY TABLE.
3. THE SURFACE OF THE TABLE IS COVERED WITH UV RESISTANT PLASTIC SHEET, WHICH IS EASY FOR MAINTENANCE, WEAR RESISTANT AND POLLUTION PROOF.
4. AG-02 EQUIPS WITH THE AUXILIARY TABLE FOR OBSTETRIC AND OTHER PURPOSE.
5. AG-02 CAN MAKE 360 DEGREE ROTATION AND CAN BE TURNED TO THE POSITION AS YOU WANT WHILE RELEASE THE TURNING LOCKER THEN LOCK IT AGAIN.



SPECIFICATIONS

1. LENGTH OF THE TABLETOP: 190 CM
2. WIDTH OF THE TABLETOP: 60 CM
3. TABLETOP ELEVATION: 65 - 95 CM
4. TRENDLENBURG AND REVERSE: 25 / 25 DEGREE
5. BACK SECTION UP/DOWN: + 60 / 0 DEGREE
6. SEAT SECTION UP/DOWN: + 25 / 0 DEGREE
7. TABLE ROTATION: 360 DEGREE

SPECIFICATIONS

- A. SHOULDER SUPPORT.
- B. HAND GRIP.
- C. PILLOW.
- D. BACK CONTROL LEVER.
- E. TREND. CONTROL HANDLE.
- F. HEIGHT ADJUSTMENT LEVER.
- G. MANUAL BRAKE.
- H. ARM BOARD.
- I. LEG CRUTCH.
- J. AUXILIARY TABLE.
- K & L: CONTAMINATION BUCKET SET.



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com

