

ESG-300 SERIES EXAMINATION UROLOGICAL AND GYNECOLOGICAL/OBSTETRIC TABLE

FEATURES:

1. FULLY AUTOMATIC ELECTRIC HYDRAULIC CONTROL SYSTEM.
2. THIS IS A MULTI-PURPOSE TABLE.
WITH ITS STANDARD ACCESSORIES, THE TABLE IS AN IDEAL EXAMINATION TABLE. WITH SOME OPTIONAL ACCESSORIES, IT CAN BE APPLIED AS A GYNECOLOGY, OBSTETRIC, PROCTOLOGY, OR URODYNAMICS TABLE.
3. **ESG-300** SERIES WITH FOUR DIFFERENT TYPE FOR YOUR CHOICE.
YOU CAN SELECT THE MOST SUITABLE ONE FOR YOUR NEED.
4. THE SURFACE OF THE TABLE IS COVERED WITH UV RESISTANT PLASTIC SHEET, WHICH IS EASY FOR MAINTENANCE, WEAR RESISTANT AND POLLUTION PROOF.
5. THE BASE IS MADE OF CAST IRON COVERED WITH HIGH QUALITY ABS PLASTIC, WHICH IS GOOD LOOKING, WATER PROOF AND EASY FOR MAINTENANCE.
6. **ESG-300** CAN OPTIONED THE AUXILIARY TABLE FOR OBSTETRIC AND OPERATION PURPOSE



FUNCTION	ESG-301	ESG-302	ESG-303	ESG-304
ELEVATING COLUMN	YES	YES	YES	YES
BACKREST ARTICULATES	YES	YES	YES	YES
SEAT INCLINES	NO	YES	YES	YES
AUTO RAISE/LOWER COLUMN	NO	NO	YES	YES
PRE-SETTING HEIGHT ADJUSTMENT (H/M/L)	NO	NO	YES	YES
AUTO TABLE AT MAXIMUM HEIGHT AND CHAIR AT MINIMUM HEIGHT	NO	NO	NO	YES


SPECIFICATIONS:

1. TABLE TOP: 550 X 1200 MM (WITH AUXILIARY TABLE): 550 X 1750 MM
2. ELEVATION: 600 - 850 MM
3. SEAT SECTION: 0 - +20 (AG-301 IMMOVABLE)
4. BACK SECTION: -15 - + 60
5. POWER SUPPLY AC 110V/220V, 240V; 50/60HZ (CAN BE FOLLOWED BY ORDER)



STANDARD ACCESSORIES:

1. SCREEN FRAME 1 PC
2. ARM BOARDS 1 PC
3. SLOPING PAN WITH FILTER 1 PC
4. LEG CRUTCHES 2 SETS
5. FORCEPS CUP 1 PC
6. FOOT CONTROLLER 1 SET
7. SEAT COVER (WATER PROOF) 1 PC
8. BED MATTRESS 1 SET



OPTIONAL ACCESSORIES:

1. DELIVERY TRAY 1 SET
2. WAIST SIDE HANDLES 2 SETS
3. AUXILIARY TABLE 1 SET
4. CONTAMINATION BASIN AND BUCKET 1 SET
5. TWIN LEG HOLDER 2 SETS
6. SPOT LIGHT (BUILT-IN) 1 SET
7. PROCTOLOGY BASIN 1 SET
8. URINE COLLECTION BASIN 1SET



ESSE 3 Via Garibaldi 30
14002 Castelnuovo D.B.(AT) - I
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.4ci.com

