



# ESG-800

**ELECTRIC MULTI PURPOSE GYNECOLOGICAL  
DELIVERY/OBSTETRIC & OPERATING TABLE**





**PATENT DESIGN**  
**EASY FOR OPERATION**  
**MULTI-PURPOSE OPERATION**

# **ESG-800**

## **ELECTRIC MULTI PURPOSE GYNECOLOGICAL DELIVERY/OBSTETRIC & OPERATING TABLE**

1. THREE-IN-ONE MULTI FUNCTION TABLE AS OBSTETRIC TABLE, GYNECOLOGICAL OPERATION TABLE, & EXAMINATION TABLE
2. ELECTRIC HYDRAULIC CONTROL FOR ELEVATION, BRAKE, BACK CHANGE AND TRENDELENBURG
3. SITTING AND LYING POSITION IS SUITABLE FOR DELIVERY. IT CAN BE AN OPERATION TABLE IMMEDIATELY WITHOUT MOVING THE PATIENT IN CASE OF NECESSARY
4. THE MAIN STRUCTURE AND ACCESSORIES ARE MADE OF STAINLESS STEEL 304 WHICH IS DURABLE AND EASY TO BE CLEANED & MAINTAINED
5. THE AUXILIARY TABLE IS CONCEALED UNDER THE TABLETOP, AND IT IS EASY TO EXTEND AND ELEVATE
6. THE TABLETOP IS AVAILABLE FOR VAGINA DOING ULTRASONIC SCANNING BY MOVING AWAY THE SLOPING PAN
7. THE FOOT CONTROLLER IS WATER PROOF AND 24V CONTROL



**TRENDELEMBURG**



**AUXILIARY TABLE APPLICATION**



**SITTING POSITION**



**LYING POSITION**




**SPECIFICATIONS:**

1. TABLETOP DIMENSION: 60X130 CM
2. TABLETOP WITH AUXILIARY TABLE 60x190 cm
3. ELEVATION: 66 - 96CM
3. TRENDELENBURG: +25°/-10°
5. BACK UP/DOWN +70°/-0°
6. POWER SUPPLY 110/220V 50/60Hz
7. OPTION : UPS WITH RECHARGEABLE BATTERIES
8. BRAKE: MECHANICAL FOOT-PEDAL

**STANDARD ACCESSORIES:**

1. PILLOW 1 PC
2. ARM BOARDS 2 PCS
3. MATTRESS 1 PC
4. WATER PROOF SEAT COVER 1 PC
5. WAIST SIDE HANDLES 2 PCS
6. LEG CRUTCHES 2 SETS
8. SLOPING PAN (WITH FILTER) 1 PC
9. FOOT CONTROLLER 1 SET
10. POWER CORD 1PC



**SCREEN FRAME**



**HEAD SIDE HANDLE**



**HEAVY TYPE LEG CRUTCHES**

## ESG 300 serie UROLOGICAL/OBSTETRIC DELIVERY TABLES

### FEATURES:

1. FULLY AUTOMATIC ELECTRIC HYDRAULIC CONTROL SYSTEM.
  2. THIS IS A MULTI-PURPOSE TABLE.
- WITH ITS STANDARD ACCESSORIES, THE TABLE IS AN IDEAL EXAMINATION TABLE. WITH SOME OPTIONAL ACCESSORIES, IT CAN BE APPLIED AS A GYNECOLOGY, OBSTETRIC, PROCTOLOGY, OR URODYNAMICS TABLE.
3. **ESG-300** SERIES WITH FOUR DIFFERENT TYPE FOR YOUR CHOICE.
- YOU CAN SELECT THE MOST SUITABLE ONE FOR YOUR NEED.
4. THE SURFACE OF THE TABLE IS COVERED WITH UV RESISTANT PLASTIC SHEET, WHICH IS EASY FOR MAINTENANCE, WEAR RESISTANT AND POLLUTION PROOF.
  5. THE BASE IS MADE OF CAST IRON COVERED WITH HIGH QUALITY ABS PLASTIC, WHICH IS GOOD LOOKING, WATER PROOF AND EASY FOR MAINTENANCE.
  6. **ESG-300** CAN OPTIONED THE AUXILIARY TABLE FOR OBSTETRIC AND OPERATION PURPOSE



FUNCTION	ESG-301	ESG-302	ESG-303	ESG-304
ELEVATING COLUMN	YES	YES	YES	YES
BACKREST ARTICULATES	YES	YES	YES	YES
SEAT INCLINES	NO	YES	YES	YES
AUTO RAISE/LOWER COLUMN	NO	NO	YES	YES
PRE-SETTING HEIGHT ADJUSTMENT (H/M/L)	NO	NO	YES	YES
AUTO TABLE AT MAXIMUM HEIGHT AND CHAIR AT MINIMUM HEIGHT	NO	NO	NO	YES

### STANDARD ACCESSORIES:

1. SCREEN FRAME 1 PC
2. ARM BOARDS 1 PC
3. SLOPING PAN WITH FILTER 1 PC
4. LEG CRUTCHES 2 SETS
5. FORCEPS CUP 1 PC
6. FOOT CONTROLLER 1 SET
7. SEAT COVER (WATER PROOF) 1 PC
8. BED MATTRESS 1 SET

### OPTIONAL ACCESSORIES:

1. DELIVERY TRAY 1 SET
2. WAIST SIDE HANDLES 2 SETS
3. AUXILIARY TABLE 1 SET
4. CONTAMINATION BASIN AND BUCKET 1 SET
5. TWIN LEG HOLDER 2 SETS
6. SPOT LIGHT (BUILT-IN) 1 SET
7. PROCTOLOGY BASIN 1 SET
8. URINE COLLECTION BASIN 1SET



ESSE 3 Via Garibaldi 30  
14022Castelnuovo D.B. (AT)  
tel +39 011 99 27 706  
fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web : [www.esse3-medical.com](http://www.esse3-medical.com)

