



# ESMC-340

## Medical Air Compressor





# ESMC-340

## Features

- CE marked, stable performance, safe and reliable
- Compact design and easy to transfer
- Assembled with clean control module, supplying cleaner and safer air to the patients
- High effective sound-absorbing materials, noise reduction
- Equipped with various ventilator of famous brand in the world

## Technical Data

- Output pressure: 0.4Mpa
- Output flow (continuous):  $\geq 40$ L/min
- Output flow (peak):  $\geq 75$ L/min
- Peak flow:  $\geq 180$ LPM
- Noise(Actually):  $< 52$ dB(A)



ESSE 3 Via Garibaldi 30  
14022Castelnuovo D.B. (AT)  
tel +39 011 99 27 706  
fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web : [www.esse3-medical.com](http://www.esse3-medical.com)

