



EGD8! \$&

Portable dental unit



DENTISTRY



EGD8! \$&

Portable dental unit



Cod. **ESPD-302**

• Technical Specification

Lightest and easy to carry multi-functional dental unit	
Low malfunction rate, low noise	
Has all functions of an integral dental unit and can be absolutely suitable for The Dental Clinic, Dental Treatment Center and the Dental Disease Prevention & Treatment Institution which can offer the door-to-door	
Flux	60L/min
Noise	40-50dB
Power supply	110V-240V 50/60Hz
Power	580W
Net weight	16kg
Dimension	350×220×450mm (LxDxH)

• Built-in Accessories

3-way syringe
4-hole(two hole) low speed handpiece hose
Four hole(two hole) high speed handpiece hose
Saliva ejector
Oil-free compressor
Clean water container 1set
Sewage water container 1set
Air container 1set

• Optional Accessories

Ultrasonic scaler
LED light curing unit
Optical handpiece
Micro-motor



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

