



ESS-TS Line

Digital Dry Incubators
with pid LCD control



STERILISATION



ESS-TS Line

Digital Dry Incubators with pid LCD control



Standard double wall fabrication.

Inner chamber made of richly anodized Stainless Steel, or highly polished stainless steel sheet.

Exterior fabricated out of thick mild sheet duly finished in staved enamel off-white paint or pastel mat finish colour combinations.

The gap between the walls is filled with-special grade glass wool for heat insulation.

Inner chamber is ribbed or placing shelves at convenient levels. Provided with two or three removable shelves.

Special elegant door closing device has been provided.

Nichrome wire heating elements provided on three sides attain quick and uniform heating in range of

50 C to 2500k \pm 10C controlled by capillary type

thermostat-Shaped thermometer is built-in type. Control panel I s provided with selector switch of high or low rates

of power, thermostat control knob and indicators for mains & thermostat. Supplied with cord and plug. The equipment

is suitable to operate on 220V AC 50 Hz single phase.

Technical Features

Cod.	Capacity	Chamber Size (HxWxD)	N. Shelves
ESS-TS28	28 L	300 x 300 x 300 mm (12"x12"x12")	2
ESS-TS45	45 L	350 x 350 x 350 mm (14"x14"x14")	2
ESS-TS65	65 L	400 x 400 x 400 mm (16"x16"x16")	2
ESS-TS95	95 L	450 x 450 x 450 mm (18"x18"x18")	2
ESS-TS125	125 L	600 x 450 x 450 mm (24"x18"x18")	2
ESS-TS224	224 L	600 x 600 x 600 mm (24"x 24"x 24")	2
ESS-TS252	252 L	900 x 600 x 450 mm (36"x 24"x 18")	2
ESS-TS336	336 L	900 x 600 x 600 mm (36"x 24"x 24")	2
ESS-TS480	480 L	900 x 900 x 900 mm (36"x 36"x 36")	3



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

