

EST-600 DELUXE MULTI-FUNCTION HYDRAULIC OPERATION TABLE



**PATENT DESIGN
EASY FOR OPERATION
MULTI-PURPOSE OPERATION**

EST-600 DELUXE MULTI-FUNCTION HYDRAULIC OPERATION TABLE

FEATURES:

1. ANTI SHOCK ABS RECOVERED STAINLESS STEEL OPERATING TABLE , FOR MULTI-PURPOSE OPERATION
2. PATENT DESIGN OF ELEVATION PILLAR OFFERS VERY STABLE TABLETOP AND IS EASY FOR OPERATION
3. HEAD AND LEG SECTION ARE INTERCHANGEABLE WHICH ALLOW FULLY C-ARM X-RAY SCANNING
4. SPECIAL DESIGN OF Y-SHAPE BASE IS GOOD FOR DOCTOR'S TREATMENT
5. SPECIAL DESIGN OF TABLETOP FRAME WITH HEAD/LEG INTERCHANGEABLE FUNCTION CAN OFFER 99% X-RAY TRANSMISSION.
6. THE TABLETOP IS MADE BY BAKELITE SHEET WHICH IS SUITABLE FOR RADIOGRAPHY AND IMAGE INTENSIFICATION PROCEDURES.
7. IT IS HARMONIOUS TO CONNECT WITH "EST-160 EXTENSION TRACTOR
8. GAS SPRING CONTROLLED BACK AND LEG SECTION PROVIDE STABLE AND SMOOTH MANIPULATION



CONNECTED WITH "EST-160 EXTENSION TRACTOR " FOR ORTHOPAEDIC TREATMENT



TRENDELENBURG POSITION

TILT POSITION

EST-600 DELUXE MULTI-FUNCTION HYDRAULIC OPERATION TABLE



KIDNEY POSITIONS



CHAIR POSITION



HEAD AND LEG SECTION INTERCHANGE FUNCTION



RADIOLUCENT TABLETOP



TREND. & TILT OPERATING

EST-600 DELUXE MULTI-FUNCTION HYDRAULIC OPERATION TABLE



1. HEAD SECTION OPERATING HANDLE.
2. BACK SECTION OPERATING HANDLE.
3. TREND./TILT OPERATING HANDLE.
4. LEG SECTION OPERATING HANDLE.
5. BRAKE PEDAL.
6. ELEVATION PEDAL.



OPTION : LEGS BECOME TWO PARTS

OPTION : SQUARE COLUMN
BECOME STAINLESS STEEL

SPECIFICATIONS:

1. TABLETOP DIMENSION: 55 X 203 CM
2. ELEVATION: 70 - 100CM (BY HYDRAULIC FOOT-PEDAL)
3. TRENDLENBURG: +30°/-30° (MECHANICAL)
4. TILT: +25°/-25° (MECHANICAL)
5. HEAD SECTION: +40°/-30° (MECHANICAL)
6. BACK SECTION: +70°/-30° (GAS SPRING)
7. LEG SECTION: +15°/-100° (GAS SPRING)
8. BRAKE: MECHANICAL FOOT-PEDAL

OPTIONS:

1. LEG PARTS BECOME TWO PARTS.
2. SQUARE COLUMN AND BASE BECOME STAINLESS STEEL.

STANDARD ACCESSORIES:

1. SCREEN FRAME.....1 SET
2. HEAD SECTION.....1 SET
3. ARM BOARD.....2 PCS
4. LEG SECTION.....1 SET
5. MATTRESS.....1 SET
6. OPERATING HANDLE.....1 PC



ESSE 3 Via Garibaldi 30
14002 Castelnuovo D.B.(AT) – I
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.esse3-medical.com

