

EST-900E DELUXE MULTI-FUNCTION ELECTRIC OPERATION TABLE



**PATENT DESIGN
EASY FOR OPERATION
FULLY ELECTRIC
MULTI-PURPOSE OPERATION
TABLE TOP SLIDING FUNCTION**

EST-900E DELUXE MULTI-FUNCTION ELECTRIC OPERATION TABLE TABLE TOP SLIDING FUNCTION

FEATURES:

1. EST-900E SERIES IS FULLY ELECTRIC MULTI PURPOSE OPERATING TABLE THE ELEVATION, TEND., TILT, BACK ARE OPERATED BY THE ELECTRIC HYDRAULIC SYSTEM.
2. THE SPECIAL DESIGN OF THE TABLE TOP CAN BE USED IN NORMAL OR INVERSE MODE, WHICH ALLOWS THE C-ARM SCAN WITHOUT BLOCKED.
3. THE EST-900E SERIES IS A MULTI-PURPOSE OPERATING TABLE, AND IT IS ALSO SPECIAL FOR THE CATH-LAB PHOTOGRAPH AND SPINE OPERATING.
4. WITH OPTIONAL ORTHOPAEDIC EXTENSION TRACTOR, THE EST-900E SERIES IS SUITABLE FOR THE USE OF ORTHOPAEDIC PROCEDURES.
5. EST-900E SERIES IS PROVIDED OF TABLE TOP SLIDING FUNCTION (30 cm)
6. COLUMN CAN OFFER 180° ROTATION
7. WIRED REMOTE CONTROL BY WATERPROOF KEYBOARD



DIRECTION CAN CHANGE OF 180°

AVAILABLE VERSIONS

"T" TYPE: WITH EXTENSION TRACTOR DEVICE FOR ORTHOPEDIC OPERATION AND FLUOROSCOPY

"R" TYPE: WITH AUTO-REVERT FUNCTION BY ONE-TOUCH FUNCTION KEY.

WITH AUTO-BRAKE FUNCTION. THE TABLE WILL NOT FUNCTION, IF IT IS NOT BRAKED.

"B" TYPE: WITH POWER AUTO RECHARGE SYSTEM FOR KEEPING THE TABLE IN POWER CONDITION

COMPOSITIONS

"EB" WITH BATTERY

"EN" WITH SPECIAL TOP ELEVATION (54-102 CM) FOR NEUROSURGERY USE

"ER" WITH AUTO-REVERSE FUNCTION

"ET" WITH EXTENSION TRACTOR DEVICE

"ERB" WITH AUTO-REVERSE FUNCTION + BATTERY

"ETB" WITH EXTENSION TRACTOR DEVICE + BATTERY

"ERTB" WITH AUTO-REVERSE FUNCTION + EXTENSION TRACTION DEVICE +BATTERY



**EST-900E DELUXE MULTI-FUNCTION
ELECTRIC OPERATION TABLE
TABLE TOP SLIDING FUNCTION**



EST-900E WITH ORTHOPEDIC EXTENSION TRACTOR AND CART



NECK FEMUR TREATMENT POSITION



TIBIA AND FIBULA TREATMENT POSITION





FEATURES

- TABLETOP DIMENSIONS 50 X 197 CM
- TABLETOP ELEVATION 74 – 110 CM (VERSION “N” : 54 - 102 cm)
- BACK SECTION UP/DOWN +90° / -40°
- LATERAL TILT RIGHT/LEFT $\pm 30^\circ$ (VERSION “N” : $\pm 22^\circ$)
- TRENDLENBURG/ANTI-TRENDELEMBURG $\pm 30^\circ$ (VERSION “N” : $\pm 22^\circ$)
- TABLE TOP ROTATION 180°
- HEAD SECTION UP/DOWN +30° / -90°
- LEG SECTION UP/DOWN +15° / -90°
- POWER SUPPLY 110V/230VAC
- POWER CONSUMPTION 500VA
- CONTROL VOLTAGE 24 VDC / 5VDC
- LOADING CAPACITY 250 kg

STANDARD ACCESSORIES:

- SCREEN FRAME \times 1 SET
- HEAD BOARD \times 1 PC
- ARM BOARD \times 1 PAIR
- FOOT BOARD \times 1 PAIR
- WASHABLE AND DEINFECTABLE SET OF MATTRESSES \times 1
- REMOTE CONTROL BOX \times 1SET



ESSE 3 Via Garibaldi 30
 14022Castelnuovo D.B. (AT)
 tel +39 011 99 27 706
 fax +39 011 99 27 506
 e-mail esse3@chierinet.it
 web : www.esse3-medical.com

