



# EST-95

*Universal Hydraulic  
Operating Table*



OPERATING THEATRE



# EST-95

## Universal Hydraulic Operating Table

### FEATURES:

- Special design for multi-purpose operation.
- Tabletop elevation is controlled by hydraulic system.
- Tabletop can be 360° rotated.
- The backrest is controlled by gas spring.
- The side rail is made of stainless steel 304.
- Accessories and clamps are made of stainless steel 304.
- The table base is covered by stainless steel 304.
- Built-in kidney elevation function, " ^ " shape.
- Brake function by footswitch.

### TECHNICAL SPECIFICATION:

- Tabletop length: 195cm
- Tabletop width: 50cm
- Tabletop elevation: 74 - 103cm
- Lateral tilt: +27° / -27°
- Trendelenburg / Reverse Trendelenburg: +30° / -22°
- Backrest vertical tilt: +60° / -20°
- Headrest vertical tilt: +15° / -90°
- Leg sections vertical tilt: +15° / -90°
- Seat section vertical tilt: +20° / -20°

### STANDARD ACCESSORIES:

- Anaesthesia frame: 1 set
- Shoulder rests: 1 pair
- Kidney bridge: 1 pair
- Arm rests: 1 pair
- Knee crutches: 1 pair
- Mattress: 1 set

### OPTIONAL ACCESSORIES:

- Fixation straps: pc
- Hand grips: pair
- U-shape headrest: set
- Raised 'UX' gHV'Y armrest: set
- X-ray radiograph sheet: set
- Contamination bucket: set



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

