

Multi-Purpose Operating Table (Side controlled)



Characteristic

- * Multi-Purpose operation table
- * Hydraulic system pump controlled
- * Tabletop can be shifted
- * Tabletop made of high strength glass, X-Ray examination can be carried out on
- * Built in kidney operation function
- * Tabletop elevation, back section and leg section are controlled by gas spring cylinder system
- * Table top shift, trendelenburg / revers, lateral tilt controlled manually
- * Side rail is made of stainless still 304

Specification

Tabletop length: 2100 mm

Tabletop width: 500 mm

Tabletop height adjusting range: 800 -1040mm

Tabletop shifting range: 400 mm

Lateral tilt range: $\geq 15^\circ$

Tredelenburg: $\geq 20^\circ$

Reveres adjustment: $\geq 15^\circ$

Head Section: Up $\geq 45^\circ$ Down $\geq 90^\circ$

Back section adjustment range: UP $\geq 45^\circ$

Leg section adjustment range: $\geq 90^\circ$ detachable

Kidney elevation range: ≥ 400 mm

Standard Accessories

Body support	1set
Arm rest	1set
Anesthesia screen	1set
Knee crutch	1set
Shoulder support	1set

EST - 160 ORTHOPEDIC EXTENSION TRACTOR

1. UNIVERSAL CONNECTION WITH ANY BRAND OPERATION TABLE.
2. INFECTION RESISTANT STAINLESS STEEL FINISHED, EASY TO BE CLEANED EVEN AGAINST THE GYPSUM APPLYING.
3. SPECIAL DESIGN FOR RADIOGRAPHY AND IMAGE INTENSIFICATION PROCEDURES.
4. AVAILABLE FOR LATERAL FEMUR, NECK FEMUR OR TIBIA & FIBULA TREATMENT POSITION.
5. THE PELVIC REST CAN BE ROTATED 360 AND IS MADE BY BAKELITE SHEET WHICH IS SUITABLE FOR RADIOGRAPHY AND IMAGE INTENSIFICATION PROCEDURES.
6. HEAVY 6 CASTERS ARE MADE BY STAINLESS STEEL, AND EACH CASTER WITH SEPARATE BRAKE DESIGN



Orthopedic Tractor



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail sse3@chierinet.it
web : www.esse3-medical.com

