



Multi-Purpose Operating Table

(Hydraulic and cylinder system)



Characteristic

- * Multi-Purpose operation table
- * Hydraulic system pump controlled
- * Tabletop made of high pressure laminate table board, X-Ray examination and fluoroscopy can be carried out on
- * Built in kidney operation function
- * Tabletop elevation, trendelenburg/reveres, lateral tilt, back section and leg section are controlled by hydraulic pump and gas spring cylinder system
- * Side rail is made of stainless still 304
- * Accessories made of stainless still 304

Specification

- * Tabletop length: 480 / 520 mm
- * Tabletop width: 1990 / 2000 mm
- * Tabletop height adjusting range: 690mm-910mm
- * Lateral tilt range: $\geq 15^\circ$
- * Tredelenburg: $\geq 20^\circ$
- * Reverses adjustment: $\geq 15^\circ$
- * Head Section: Up $\geq 35^\circ$ Down $\geq 40^\circ$ detachable
- * Back section adjustment range: UP $\geq 75^\circ$ Down $\geq 15^\circ$
- * Leg section adjustment range: Down $\geq 90^\circ$ detachable (bifurcation angle $\geq 180^\circ$)
- * Kidney elevation range: ≥ 80 mm

Standard Accessories:

Cassette rack	1set
Body support	1set
Single-sheet arm rest	1set
Handle of kidney bridge	1set
Anesthesia screen	1set
Knee crutch	1set
Body strap	1set

Optional Accessories:

Double-sheet arm rest
Versatile head rest
Tractor for orthopedics
Head rest



ESSE 3 Via Garibaldi 30
14022Castelnuovo D.B. (AT)
tel +39 011 99 27 706
fax +39 011 99 27 506
e-mail esse3@chierinet.it
web : www.esse3-medical.com

