



# FD306

Pocket foetal doppler



INTENSIVE CARES



# 8' \$\*

## Pocket foetal doppler



- LCD backlight display with:
  - FHR data
  - Probe working frequency
  - Battery charge
  - Doppler signal
  - Alarm indicator
- User friendly interface for easy operation
- Amplified sound and recorded function
- Clear sound with a high quality loudspeaker
- 2.0MHz standard probe

Cod. F8' \$\*

### Technical Features

Ultrasonic frequency	2.0 MHz(with the standard probe)
Ultrasonic power intensity	≤10mW/cm <sup>2</sup>
Comprehensive Sensitivity	≥90dB
Measurement range	50bpm-210 bpm
Isptp	≤0.1Mpa
Doppler Frequency	0-3KHz
Power Consumption	≤1.5W
Automatic power off	after 3 minutes
Working Mode	Continuous ultrasound Doppler
Working temperature	5-40°C
Humidity	<80%
Atmospheric pressure	86-106kPa
Power Supply	3 x 1,5V battery
Optional available probes	3.0MHz, 4.0MHz, 5.0MHz and 8.0MHz



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)



CE