



ITALY & GERMANY

FXD-150+P

Flexible use
150watt Xenon light source with pump



FXD-150+P

FXD-150+P offers the professional user the most advanced pure day light quality sources(5700°K color temperature). It contains an improved high intensity and low heat conductive xenon lamp.

Lamp:

The output intensity is controlled stepless via a motor driven mechanical diaphragm. It ensures a constant color temperature. The output intensity is set by push buttons.

Connection:

An adapter for ESSE 3 & Wolf. rigid endoscopes application. Flexible fiberscope connection which is compatible with Olympus.

Pump:

With integrated pump.

Quality:

The light source fulfills all essential safety demands. The rugged desk housing and the flat, waterproof touch control panel withstands all daily used requirements.

Brightness indication:

With the LED indicator for the brightness, makes brightness more accurate.

Lamp Span indication:

With LED lamp span indicator, the Xenon Lamp life can be easily recognized.

TECHNICAL DATA

Model: FXD-150+P

Power supply: AC(110-220V)±10%,50-60HZ

Dimensions: 350x120x370

Weight: 12KG

Power consumption: 200WATT

Adapter : ESSE 3& Olympus

With integrated pump

Life span: 500 hours

Brightness: 570000x

Lamp: 150Watt Xenon light lamp

Manufactured and tested: ISO9001/2000, CE

ESSE 3 Medical SRL
Via Garibaldi 30
14022 castelnuovo DB
Tel +39 011 99 27 706
Fax +39 011 99 27 506
www.4ci.com