



# Cylinder Data Chart

Medical Gas/ Medical Gas Mixture	Colour code	Cylinder code	C	AD*	DD*	PD	RD	SD*	D*	E	AF	DF	F*	LF	VF	AV	HX	T4*	T7*	G*	AK	J	L
		Approx dimensions mm in	430 x 89 17 x 3.5	480 x 100 19 x 4	480 x 100 19 x 4	455 x 100 18 x 4	480 x 100 19 x 4	430 x 100 17 x 4	535 x 102 21 x 4	865 x 102 34 x 4	670 x 175 26 x 7	690 x 175 27 x 7	930 x 140 37 x 5.5	930 x 140 37 x 5.5	930 x 140 37 x 5.5	795 x 181 31 x 7	940 x 140 37 x 5.5	1205 x 195 47 x 8	1205 x 195 47 x 8	1320 x 178 52 x 7	1640 x 229 65 x 9	1520 x 229 60 x 9	1640 x 230 65 x 9
		Approximate weight (empty) kg lb	2.0 4.5	3.7 8.0	4.1 9.0	4.8 10.0	4.1 9.0	3.1 7.0	3.4 7.5	5.4 12.0	9.9 22.0	10.0 22.0	14.5 32.0	14.5 32.0	14.5 32.0	15.5 34.1	17.4 38.3	32.5 72.0	32.5 72.0	34.5 76.0	51.0 112.0	68.9 152.0	51.0 178.0
<b>OXYGEN</b>		Contents (litres) Valve type	170 Pin-index	400 Star valve 6mm tubing	460 Star valve 6mm tubing	300 Bullnose 5/8" BSP (f)	460 Star valve Schrader/ 6mm tubing	300 Pin-index	340 Pin-index	680 Pin-index	1360 Bullnose 5/8" BSP (f)	1360 Star valve Bayonet	1360 Bullnose 5/8" BSP (f)				2300 Star valve Schrader			3400 Bullnose 5/8" BSP (f)		6800 Pin-index (side spindle)	
<b>NITROUS OXIDE</b>		Contents (litres) Valve type	450 Pin-index						900 Pin-index	1800 Pin-index			3600 Handwheel 1 1/16" x 20tpi (m)							9000 Handwheel 1 1/16" x 20tpi (m)		18000 Handwheel 1 1/16" x 20tpi (m)	
<b>ENTONOX (50% N<sub>2</sub>O/50% O<sub>2</sub>)</b>		Contents (litres) Valve type							500 Pin-index				2000 Pin-index								5000 Pin-index		
<b>AIR</b>		Contents (litres) Valve type							640 Pin-index				1280 Bullnose 5/8" BSP (f)					4000 Star valve Schrader	4000 Star valve Schrader	3200 Bullnose 5/8" BSP (f)		6400 Pin-index (side spindle)	
<b>CARBON DIOXIDE</b>		Contents (litres) Valve type	450 Pin-index						1800 Pin-index				3600 Handwheel 0.860" x 14tpi (m)	3600 Handwheel 0.860" x 14tpi (m)									
<b>OXYGEN/CARBON DIOXIDE MIXTURE (95% O<sub>2</sub>/5% CO<sub>2</sub>)</b>		Contents (litres) Valve type											1360 Bullnose 5/8" BSP (f)							3400 Bullnose 5/8" BSP (f)		6800 Bullnose 5/8" BSP (f)	
<b>HELIUM</b>		Contents (litres) Valve type							300 Pin-index				1200 Bullnose 5/8" BSP (f)										
<b>HELIUM/OXYGEN MIXTURE (79% He/21% O<sub>2</sub>)</b>		Contents (litres) Valve type											1200 Bullnose 5/8" BSP (f)										
<b>LUNG FUNCTION MIXTURES TYPE 1-4</b>		Contents (litres) Valve type/Outlet type															1500 Side outlet Handwheel 5/8" BSP (LH) (f)				6000 Side outlet Handwheel 5/8" BSP (LH) (f)		
<b>CARBON DIOXIDE/OXYGEN MEDICAL GAS MIXTURES</b>		Contents (litres) Valve type															1460 Side outlet Handwheel 5/8" BSP (f)						7300 Side outlet Handwheel 5/8" BSP (f)
<b>CARBON DIOXIDE/AIR MEDICAL GAS MIXTURE (5% CO<sub>2</sub>/95% AIR)</b>		Contents (litres) Valve type															1350 Side outlet Handwheel 5/8" BSP (f)						6750 Side outlet Handwheel 5/8" BSP (f)
<b>HELIUM/OXYGEN/NITROGEN MEDICAL GAS MIXTURE (56% N<sub>2</sub>/35% O<sub>2</sub>/9% He)</b>		Contents (litres) Valve type															1310 Side outlet Handwheel 5/8" BSP (f)						6550 Side outlet Handwheel 5/8" BSP (f)
		<b>Water capacity (litres)</b>	1.2	2.0	2.0	2.0	2.0	2.0	2.32	4.68	9.43	9.43	9.43	9.43	9.43	10.0	10.0	20.0	20.0	23.6	40.0	47.2	50.0



NOTES: \* The indicated cylinders are for specialised applications and availability is restricted.

- Most common cylinders used within the hospital.
- Specialist use cylinders used within the hospital.
- Other medical cylinders used within the hospital.

ESSE 3 Via Garibaldi 30  
14022Castelnuovo D.B. (AT)  
tel +39 011 99 27 706  
fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web : [www.esse3-medical.com](http://www.esse3-medical.com)