



High-end
video - surgery system
with a harmonious interplay
of the individual components.



ENT

The symbiosis of medical
knowledge and engineering
know-how




VIDEOSYSTEMS ENT

ENT

ENT

AUTOCLAVABLE 134° C / 273° F



Rigid ENT

SINUSCOPES

Diameter	2.7mm	2.7mm	2.7mm
Length	100mm	110mm	110mm
Direction of View	0°	0°	30°
Item-No.	521001	521002	521003


Diameter	4mm	4mm	4mm	4mm	4mm	4mm
Length	50mm	110mm	110mm	175mm	175mm	175mm
Direction of View	0°	0°	30°	0°	30°	45°
Item-No.	521006	521008	521009	521015	521016	521017

Diameter	4mm	4mm
Length	175mm	175mm
Direction of View	70°	90°
Item-No.	521018	521019

ENT

OTOSCOPES

AUTOCLAVABLE 134° C / 273° F



Rigid OTOSCOPES


OTOSCOPES

Diameter	4mm	4mm	4mm	4mm	3mm	3mm
Length	50mm	90mm	90mm	90mm	90mm	90mm
Direction of View	0°	0°	30°	70°	0°	30°
Item-No.	521020	521021	521022	521023	521024	521025

Diameter	2.7mm	2.7mm	2.7mm	2.7mm	2.7mm	2.7mm
Length	90mm	90mm	90mm	110mm	110mm	110mm
Direction of View	0°	30°	70°	0°	30°	70°
Item-No.	521026	521027	521028	521029	521030	521031

Laryngo-Pharyngoscopy

Laryngo-Pharyngoscope



Sharp image, significant depth of field

- High performance illumination
- 2 different magnifications:
Panoramic view and examination of details
- Soakable and gas sterilizable


Magnification type	
Diameter	10mm
Length	153mm
Viewing direction	90°
Working distance	15mm
Magnification	1 : 1.6
Item-No.	521060

Diameter	6mm	6mm	8mm	8mm
Length	185mm	180mm	185mm	180mm
Direction of View	70°	90°	70°	90°
Item-No.	521050	521051	521053	521054



ENT

ENT Sheath with rotating stopcocks


ENT Sheath with two rotating stopcocks





AUTOCLAVABLE 134° C / 273° F

Scope Ø	Scope Length	 Item-No.:	 Item-No.: 30°
2,7 mm	175 mm	522000	522001
2,7 mm	187,5 mm	522002	522003
4,0 mm	175 mm	522004	522005
5,0 mm	175 mm	522006	522007

ENT Sheath with one rotating stopcock



AUTOCLAVABLE 134° C / 273° F

Scope Ø	Scope Length	 Item-No.: 0°	 Item-No.: 30°
2,7 mm	175 mm	522012	522013
2,7 mm	187,5 mm	522014	522015
4,0 mm	175 mm	522016	522017
5,0 mm	175 mm	522018	522019

ENT

Trocars and Obturators for ENT Sheaths

Trocar pyramidal



AUTOCLAVABLE 134° C / 273° F

Scope Ø	Scope Length	Item-No.:	
2,7 mm	175 mm	522020	
2,7 mm	187,5 mm	522021	
4,0 mm	175 mm	522022	
5,0 mm	175 mm	522023	

Obturator blunt



AUTOCLAVABLE 134° C / 273° F

Scope Ø	Scope Length	Item-No.:	
2,7 mm	175 mm	522020-b	
2,7 mm	187,5 mm	522021-b	
4,0 mm	175 mm	522022-b	
5,0 mm	175 mm	522023-b	

Obturator blunt



AUTOCLAVABLE 134° C / 273° F

Scope Ø	Scope Length	Item-No.:	
2,7 mm	175 mm	522024	
2,7 mm	187,5 mm	522025	
4,0 mm	175 mm	522026	
5,0 mm	175 mm	522027	



ENT

Sheaths with rotating stopcocks for small joints, Highflow - Stainless Steel - Tubes drawn in

ENT Sheath with two rotating stopcocks





AUTOCLAVABLE 134° C / 273° F

Scope Ø	Scope Length	Outer Ø	 Item-No.: universal	 Item-No.: 30°
2,7 mm	175 mm	4,5 mm	522028	522029
2,7 mm	187,5 mm	4,5 mm	522030	522031
4,0 mm	175 mm	6,0 mm	522032	522033

ENT Sheath with one rotating stopcock



AUTOCLAVABLE 134° C / 273° F

Scope Ø	Scope Length	Outer Ø	 Item-No.: universal	 Item-No.: 30°
2,7 mm	175 mm	4,5 mm	522034	522035
2,7 mm	187,5 mm	4,5 mm	522036	522037
4,0 mm	175 mm	6,0 mm	522038	522039

ENT

Trocars and Obturators for ENT Sheaths, Highflow- Stainless Steel-for Tubes drawn in

Trocar pyramidal



Scope Ø	Scope Length	Item-No.:
2,7 mm	175 mm	522040
2,7 mm	187,5 mm	522041
4,0 mm	175 mm	522042

Obturator blunt

AUTOCLAVABLE 134° C / 273° F



Scope Ø	Scope Length	Item-No.:
2,7 mm	175 mm	522043
2,7 mm	187,5 mm	522044
4,0 mm	175 mm	522045

Obturator conical blunt

AUTOCLAVABLE 134° C / 273° F



Scope Ø	Scope Length	Item-No.:
2,7 mm	175 mm	522046
2,7 mm	187,5 mm	522047
4,0 mm	175 mm	522048

Endoscope Handle, round



Item-No.: 524001

Gumholeaspirator, with short edge



Item-No.: 524002

Suction and Irrigation Tube



Item-No.: 524003


Suction Canula, short, smooth



Item-No.: 524004

Antrum- and Sinuscopy Trocar


Ø 5.0 mm, incl. pyramid-shaped trocar



Item-No.: 524005

Biopsy Punch, straight, vertical

3.0 mm double mobile



Item-No.: 524006

Nasal Cutting Forceps, Weil-Blakesley,

sharp, without neck, straight, 190 mm



Item-No.: 524007

Nasal Cutting Forceps, Weil-Blakesley,

sharp, without neck, 45°, 190 mm



Item-No.: 524008

Nasal Scissors, straight




Item-No.: 524009

Nasal Scissors, right curved




Item-No.: 524010

Sickle Knife, 200 mm




Item-No.: 524011

Sickle Knife, 200 mm, long curved




Item-No.: 524012

Spoon, sharp, round form, 200 mm



Item-No.: 524013

Spoon, sharp, round form, 200 mm




Item-No.: 524014

Spoon, sharp, oval form, 200 mm



Item-No.: 524015

Septum-Elevator, Freer
sharp/blunt, 180 mm



Item-No.: 524016



ESSE 3 Medical
Via Garibaldi 30
14022 Castelnuovo DB
ITALY
tel +39 011 99 27 706
fax +39 011 99 27 506
www.esse3-medical.com