



MADE IN GERMANY & JAPAN

Semi flexible shaft
Distal end with atraumatic form
High resolution picture bundle
Fix Focus
Optional holder for ultrasonic or
Electro hydraulic leads, eg EMS
Telescopic bridge for Wolf

sapphire glass lenses
brilliant, crystal-clear picture
absolutely tight sealing
non-corrosive
ergonomic design
134 °C/270 °F autoclave



URETERO - RENOSCOPE
CYSTOSCOPE-FLEXIBLE


URETERO-RENOSCOPE

URETERO-RENOSCOPE





URETERO RENOSCOPES

Made entirely from non-corrosive, stainless steel.
 All mechanical connections are laser-welded. Distal and proximal lenses are made of sapphire glass and soldered in place.
 This means no joins offer brilliantly clear picture quality.
 ESCAD Optics are compatible with the shafts of other brands and they perfectly integrate with your existing equipment.
 Optics are with a range of color coded viewing angles.
 The blue ring stands for autoclave.



Compact uretero renolescope with side view



Compact uretero renolescope



Distal
6.0 Charr.

First step
7.5 Charr.












Distal
7.5 Charr.

First step
9.8 Charr.

Uretero Renoscope, Semi-flexible, Up to 15degree (Laterally)						
Distal	First step	Working length	Working channel	View direction	Angle	Item No.
Ch 6	Ch 7.5	310mm	1x4 or 2x2.4 Ch	8 degree	85 degree	916004
Ch 6	Ch 7.5	425mm	1x4 or 2x2.4 Ch	8 degree	85 degree	916005
Ch 7.5	Ch 9.8	310mm	1x5 or 2x3 Ch	8 degree	85 degree	916006
Ch 7.5	Ch 9.8	425mm	1x5 or 2x3 Ch	8 degree	85 degree	916007

URETERO-RENSCOPE



Accessories for Uretero Renoscope, semi flexible, 60cm length				
		Item No.		
		Ch 3	Ch 4	Ch 5
	Biopsy forceps	916100	916101	916102
	Foreign body biopsy forceps	916103	916104	916105
	Alligator forceps	916106	916107	916108
	Rat Tooth forceps		916109	916110
	Stone street forceps	916111	916112	916113
	Ureter Strictures scissors		916114	916115
	Three prong retriever	916116		916117
	Dormia stone catcher	916118		916119
	Single action scissors	916120		916121



CYSTOSCOPE-FLEXIBLE

Made entirely from OLYMPUS IMAGE BUNDLE.
 All mechanical connections are super combined. Distal and proximal lenses are made of super clear image bundle from Japan.
 This means no fog after certain time use.
 and provide always brilliantly clear picture quality.



Fully water-proof secures long term use.

CYSTOSCOPE-FLEXIBLE	
Item No.	ES -CY 4
Distal Diameter	4.9mm
Instrument Channel	2.0mm
Forward View	100 Degree
Depth of View	3.50mm
Bending section	Up 210 Degree , Down 100 Degree
Working length	380mm

Scope of delivery:
 case, valves, o-rings for sealing, leakage cap, cleaning brush short/long,
 cell brush, biteblock, biopsy forceps - oval cup without spike, leakage tester
 cleaning adapter, siliconoil, manual