

PRODUCT DESCRIPTION



ITALY & GERMANY

P 750II

AUTO PUMP SUCTION & IRRIGATION



[www.4ci.com](http://www.4ci.com)

## Suction & Irrigation pump P750-II

The irrigation power of the pump is automatically operated by a restricted pressure sensor.

The regulation of Irrigation and Suction is automatically controlled and restricted by integrated micro-computer sensor .

This unique design is more advanced by our original German model.

The pump is very simple to install and easy to use.

## TECHNICAL DATA

Power supply:	AC220±10% switchable
Irrigation rate:	0---1.5L/ Min
Suction rate:	2.5L / Min
Dimension :	380X310X130MM
Max pump capacity:	750mmHg
Min vacuum suction:	-75kpa
Weight:	7kg
Power consumption:	100 W
Working noise:	≤ 55DB
Class/type :	I/BF
Working atmosphere:	0C°-40 C°
Manufactured and tested according: ISO13485&ISO9000,CE	

## Delivery contains

P750II Pump

silicone tube set

**ESSE 3 Export SRL**  
Via Garibaldi 30  
14022 castelnuovo DB  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
[www.4ci.com](http://www.4ci.com)