

PRODUCT DESCRIPTION



ES-APU-10
Arthroscope Pump system



www.4ci.com

ES-APU-10

The shaver system is also a complete micro-computed controlled unit. Combination working with WKS Arthroscope shaver, the unique interfrated pressure sensing system provides instantaneous safety, precise adjustment throughout the whole operation procedure.

- Adjustable speed from 200 to 500 RPM.
Fully meet all arthroscope field operation's requirement..
- For use during interventions on small and large joints.
- Long noise during operation.
- With equipped accessories secure sterilized operation.
- Calm and smooth motion, high suction capacity.
- User-friendly design to make the operation efficient.

TECHNICAL DATA

Power supply:	230 V AC- 50 Hz
Power consumption:	90 VA,
Speed:	200-500 rpm
Pressure:	0-280 mm Hg
Pressure control:	electronical
Dim.: (WxHxD):	112x400x300mm
Protective class:	I , Typ BF, IPX0,
EN-Class:	2a,
Weight:	17 kg



ESSE 3 Export SRL

Via Garibaldi 30
14022 castelnuovo DB, Italy
tel +39 011 99 27 706
fax +39 011 99 27 506
esse3@chierinet.it
info@4ci.com
www.4ci.com