







Object : MOBILE MEDICAL CT SCANNER UNIT

		Qté	
A	<p>HITACHI Whole Body X-ray CT System SUPRIA 16 ch. / 32 slices with 5.0MHU X-ray tube) < configuration per set ></p> <p>1) Scanner gantry ---- 1 2) Patient table ---- 1 3) Operator's console ---- 1 + Console unit (1) + 24" Wide LCD monitor (1) + Keyboard (1) + Mouse (1) + Interphone box (1)</p> <p>4) Standard functions ---- 1 + DICOM 3.0 (1) + DICOM Print (1) + Predict Scan (1) + Intelli-IP (Advanced) (1) + Intelli-EC (SD-mode) (1) + CEV-CPR (1) + Auto MPR (1)</p> <p>5) Standard accessories ---- 1 + Head rest (without arm rest) (1) + Head/Chin band (1) + Fixing band (1) + Patient mat (1) + Patient belt (L)/(R) (1) + Speaker (for CT room) (1)</p> <p>6) Foot switch ---- 1 7) Instruction manual ---- 1 8) System transformer (*1) ---- 1 9) Tera Recon Post-Processing (including 6 advanced modules) Workstation : Vascular, Dental, Virtual Endoscopy/Coloscopy, CAD Lung and colon, Segmentation Analysis & Tracking, Finding Workflow – RECIST) with workstation and 24" LCD 10) Automatic Injector (single-head) 11) Laser Imager (on dry films) HighCapXp 12) UPS 80 kVA – 10 mn 13) Installation (by trained staff)</p> <p>Power requirement (main system) : 3 phase, 380V, 75kVA (*1) : provided if the power supply voltage at site is 380/400/415/440 VAC</p>	 <p>Supria 16ch/32slice OPEN & COMPACT LOW DOSE CT</p>	1
	B < OPTIONS for SUPRIA >		
	1)	Breath-Navi display	1
	2)	Light localizer	1
	3)	Arm Rest HF + Head Rest2 + Chin Rest	1
	4)	Arm Rest FF + Head Rest3	1
	5)	Spacer1 (angle 5deg.)+Spacer2 (angle 10deg.)	1
	6)	Triangle Mattress + Foot Mattress FF	1

7)	Infant Fixing Tool	1
8)	Arm Rest for intravenous drip	1
9)	DICOM MWM	1
10)	DICOM MPPS	1
11)	DICOM Query/Retrieve	1
C	<p>Electric generator 80KVA ES-HHC/S MANUAL/automatic ELECTRIC START</p> <p>The group with soundproof canopy will be housed inside the medical cabin in a space that is insulated and soundproof, easy to inspect; the rest of the available space (width 2.40 and height 2.40 mt and depth 1 mt) will be used as a storage cabinet for accessories and films</p> <p>Equipped with sound proof canopy, water jacket heater, 80KVA protection inverter, complete with rechargeable battery for manual and automatic start and automatic transfer switch (in case of conenction directly to an external electric power supply net)</p>	1
D	<p>SHELTER 40" TYPE HQ External height 2896 mm, internal 2.698 mm External length 12192 mm, internal 12032 mm External width 2438 mm, internal 2352 mm Curb weight 3,940 kg Capacity 26540 kg Capacity 76MC</p>	1
E	<p>INTERNAL INSTALLATION With</p> <ul style="list-style-type: none"> - Soundproof service section for <ul style="list-style-type: none"> ● the generator set ● the automatic start cabinet ● 80KVA UPS ● the air conditioning unit <p>Equipped with internal double accessing door and internal leaded inspection door</p> <ul style="list-style-type: none"> - Observation / office section with leaded glass 100x80 cm <p>Equipped with external accessing door and internal leaded door</p> <ul style="list-style-type: none"> - CT scanner section with full leaded pavement and walls <p>Equipped with 2 external maintenance/inspection leaded doors</p> <ul style="list-style-type: none"> - Electrical system for the 3 sections - Air conditioning of both sections - Complementary accessories (entrance door, removable ladder, retractable awning, accessing 	1
F	<p>RADIOPROTECTION MATERIAL</p> <ul style="list-style-type: none"> - Qty 10 rolls of leaqd for pavement, walsl and door of radiology section <p>set of radioprotection accessories</p> <ul style="list-style-type: none"> - 1 leaded aron - 1 pair of leaded gloves - 1 gonad protection set - 1 ovarian protection set - 1 thyroid protection set - 1 leaded window 	1
G	<p>LOW SEMI TRAILER for 40' container</p> 	1
Option		
H	<p>TRACTOR TRUCK IVECO STRALIS <i>See ANNEXE I</i></p> 	



Supria ^{16ch/32slice}
 OPEN & COMPACT LOW DOSE CT



START!

WIDE & COMPACT EASY & QUICK

ANNEXE I

If you want purchase the truck from us, we can propose

- Truck **IVECO STRALIS** , offered in 2 versions
Brochure déchargeable sur https://issuu.com/iveco1975/docs/new-iveco-stralis-brochure_b7a8da28b0c546?e=5200183/38765935
 - **IVECO AS440X46T HR ON+**(mechanical suspensions, *delivery time 5/6 months*)
 - **IVECO AS440S46T/P** (pneumatic suspensions, *delivery time 2 months*)

1st proposition **IVECO AS440X46T HR ON+**

AS HYWAY cabin

Clutch: manual 16 speeds

Axis configuration: 4x2

Front / rear disc brakes

Drive to the left

MTTG VW or GCW 18-20 tons

3800 mm pitch

Bridge: double reduction

Engine power C11 460 hp

Mechanical suspensions

Version: tractor

Equipped with

- Vehicle Kipper
- Reinforced front parabolic springs
- Lighthouse protection grille
- Front bumper made of sheet metal
- Ralla of 2 height 150 mm + plate 50 mm
- Ralla JOST
- Bridge report 4.23
- Spare wheel
- Inflating hose Tire 15 meters
- Tires 315/80 R22,5 versatility of use
- Battery 220Ah
- ABS + ASR brakes
- Cabin climate control with automatic adjustment
- Refrigerator cabin 50 lt
- 2 bunk beds with cab protection
- Leather steering wheel

Delivery time : 5/6 months



2nd proposition **IVECO AS440S46T/P**

AS HYWAY cabin

Clutch: HI TRONIX 12 speed

Axis configuration: 4x2

Front / rear disc brakes

Drive to the left

MTTGVM or GCW 18-20 tons

3800 mm pitch

Bridge: simple reduction

Engine power C11 460 hp

Pneumatic suspensions

Version: tractor

With

- Telematic box
- Driving style evaluation + driver attention support (DSE + DAS)
- Vehicle Data Interface (FMS = Fleet Meating system)
- Iveconnect (touch screen radio with browser)
- Refrigerator cabin 50 lt
- HIGH CONFORT pneumatic driver's seat in fabric (adjustable lumbar, ventilated, heated)
- Passenger seat in pneumatic adjustment fabric
- Tire Pressure Measurement Instrument (TPMS)
- Measure loads on the axes
- Spare wheel
- Ralla of 2 height 150 mm + plate 50 mm
- Ralla JOST
- 770 lt aluminum tank
- Predictive GPS DRIVE
- Cabin air suspension
- Cabin climate control with automatic adjustment
- 1 lower bed with compartment and cab protection
- Inflating hose Tire 15 meters
- Tires 315/80 R22,5 versatility of use
- Battery 220Ah
- ABS + ASR brakes
- Integrated hydraulic retarder
- Upper and lateral spoiler

Delivery tive : 2 months

