



ESCH-2000

DENTAL ULTRASONIC PIEZO SCALER MACHINE



Feature :

Multiple function

with different scaler tips , the machine is applicable to implant maintenance , broken file removing , gutta -percha cutting , crown removing , cavity preparation , shoulder preparation , calcifications removing and root planning .

Sensitive touch screen control panel , steel glass water proof

High efficient root canal irrigation.

Smart touch panel control

Self water supply

Anti-corrosion internal pipeline : Hydrogen peroxide .Sodium hypochlorite and other special chemicals are all applicable .

Equipment list

Handpiece 1

Tips 8

Torque wrench 1

600ML water tank 1

Technical parameters

Torque	10 power
Tips model	G1/G2/G4/G5/G6/ P1/P3/P4/P4D/E1
Working frequency	28kHz+3kHz
Amplitude	20-60µm

Key attributes

Industry-specific attributes

Woodpecker compatible cavitron piezo ultrasonic scaler LED detachable handpiece self water supply for scaling perio root canal endo irrigation

Other attributes

Instrument classification: Class II

Power Source: Electricity

Quality Certification: CE

Product name: Ultrasonic scaler LED

Certificate: CE ISO

Function: Perio root canal endo irrigation scaling

Feature: With different tips, can do many jobs

Handpiece: LED detachable handpiece

Display: Touch control panel

Self water system: Yes ,600ML self -watertank supply

Water pressure: 0.1bar-5bar(0.01Mpa-0.5Mpa)



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail: danielec@esse3-medical.com
web: www.esse3-medical.com

