



DATA SHEET

Model: CP-800

Ophthalmic chart projector

Can provide 30 types of clear icons with different functions, and can program icon groups per specific needs. It can produce more precise measurement results in shorter time by remote controller and fast toggling icon, applying red/green & polarized icon test sheet, binocular eyesight and binocular balance..

Easy installation and maintenance

It is quite easy to replace light bulbs, and unnecessary to adjust its focus and position

Customer Programmable

Icons can be displayed in predefined sequence by programming.

Display Single Icon

Display single icon by multi-function mask plate.

Tilted Placement

Perfect image position can be obtained by adjusting the horizontal axis of projector



Projection distance: 1.5m~ 6m

Projection magnification: 30 x (at 5m)

Projection size: 330mm(W) x270mm(H) (at 5m)

Chart: 30 different pattern

Speed of chart conversion: One chart per 0.03s

Mask 1 open, 5 horizontal lines,8 vertical lines, 21 single letters, 1 red/green

speed of mask conversion: One mask per 0.03s

Program: 2 sets programs, each program contains up to 30 steps

Light source: 12V,50W Halogen Lamp/ LED LAMP

Auto-off function: After 10 minutes idle time.

Power source: AC 220V ,50Hz or 110V, 60Hz


Power consumption: 80W

Dimensions: 300mm(L) x 230mm(W) x 240mm(H)

Net weight: 6Kg

Accessories: Remote control, Polarized metal screen, halogen lamp, polarized glasses, fuses(2), batteries(2)

Optional accessories: Floor stand, wall standard



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

