



ES-AUN Series

Bench Top Steam Sterilizer



STERILIZATION



ES-AUN Series

Bench Top Steam Sterilizer



Characteristic

- * Automatic and safety door lock system
- * LCD display, easy for operation
- * Steam compensation techniques
- LCD display. 2 programs, easy to choose, 121/134.
- Fully-Automatic feeding water, excess water and air to be drained away
- Display the failure reason on LCD screen. Easy to correct.
- Can show the reason of failure on LCD, don't need to check recording to code menu.

Safety Devices

- Machine has safety door lock system
- Over heat control: when inner temperature is over heated, power will be switch off
- Automatically and start alarm system
- Over pressure control: alarm will be announced when the machine is over pressure
- Low water level detected: alarm will be announced when it is low water

Specifications

Model	ES-AU14N	ES-AU17N	ES-AU23N	ES-AU30N	ES-AU45N
Chamber (mm) (Φ×D)	250x280mm	250x380mm	250x490mm	250x600mm	350x650mm
Chamber Capacity (liter)	14L	17L	23L	30L	45L
Overall Size(cm) (L×W×H)	48x48x42	48x53x42	48x69x42	80x48x42	90x60x55
Net Weight (KGS)	31	34.5	41	53	65
Gross Weight (KGS)	38	43	50	62	79
Packing Size(cm)	68x60x50	73x60x50	87x60x50	96x57x50	110x71x71
Power	AC220V, 110V, 126V	AC220V, 110V, 126V	AC220V, 110V, 126V	AC230V, 50/60Hz	AC220V±10%, 50Hz±10%
Heater Supply	1200W	1200W	1500W	1500W	2000W
Sterilize Temperature	121/134°C	121/134°C	121/134°C	121/134°C	121/134°C
Sterilize Time	Please refer to different programs				
Dry Procedure	Vacuum drying				
Sterilizer Construction	Stainless steel SUS 304				
Water tank					2 water tanks inside
Stainless steel instrument contain	2 trays	2 trays	2 trays	1 tray holder, 5 mesh trays	3 trays



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

