

**Cod. NB-50B**




**Cod. NB-50T**

Stereomicroscope line designed for operators that require high reliability and precision, finding space even in research laboratories. The accessories will increase the potential.

### TECHNICAL SPECIFICATION

Cod.	NB-50B	NB-50T
Binocular head inclined 45°	✓	OPTIONAL
Trinocular head inclined 45°	OPTIONAL	✓
Extra wide field eyepiece EWF10x/22mm	✓	✓
Diopter adjustment (±5mm) - Interpupillary distance (55-75mm)	✓	✓
Zoom objective 0.8x....50x	✓	✓
Zoom ratio 6.3:1	✓	✓
Working distance 115mm	✓	✓
Focusing system with coaxial coarse	✓	✓
LED Illumination incident adjustable	✓	✓
LED Illumination transmitted adjustable	✓	✓
Glass insert opal & black acrylic plate (Ø125mm)	✓	✓
Couple rubber ring	✓	✓
Dust cover	✓	✓
Network cable	✓	✓
Foto/Video Adapter	/	✓



### Configuration for application in Electronic/Mechanical:

<b>SV-6</b>	Stand for microscope
<b>NB-TH</b>	Trinocular head zoom
<b>NB-FS</b>	Sistema di messa a fuoco

The head can be NB-BH binocular.



### Accessories:

<b>NB-BH</b>	Binocular head zoom 0.8x - 5x with eyepiece.
<b>NB-TH</b>	Trinocular head zoom 0.8x - 5x with eyepiece
<b>NB-FS</b>	Sistem focus for NB/BH and NB/TH head.
<b>B2-1500</b>	Micrometer eyepiece WF 10x
<b>B2-1505</b>	Couple of eyepieces WF15x/16mm
<b>B2-1510</b>	Couple of eyepieces WF20x/12mm
<b>B2-1515</b>	Additional lens 0.5x
<b>B2-1520</b>	Additional lens 0.7x
<b>B2-1525</b>	Additional lens 2x
<b>B2-1530</b>	CCD adapter 1x with C-mount (for trinocular)
<b>B2-1535</b>	CCD adapter 0.5x with C-mount (for trinocular)
<b>B2-1540</b>	Plastic ring for eyepiece (couple)
<b>B2-1545</b>	Polarizing set
<b>B1-0170</b>	Micrometer cover glass
<b>B1-0175</b>	Bag for carrying



**ESSE3** srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

