




MINI Line

Electronic Scales



LABORATORY SCALES



MINI Line

Electronic Scales

Electronic scale with good performance, small, compact and easy to handle and to use.

Technical specification:

- LCD display with backlight.
- Auto power-off.
- External calibration.
- Piece counting.
- Stability indicator.
- Automatic tare.
- Selectable units: g,Kg,oz,lb.
- Battery level indicator.
- Power supply AC 12V/250mA.
- Operates with 4 AA batteries (not included)
- Working temperature 0 - 40 °C.
- Net weight 650 g.
- Dimensions: 180x230x50h



Cod.	Capacity	Resolution	Linearity	Pan size
MINI 500	1000 g	0,1 g	± 0,1 g	155X155 mm
MINI 2000	3000 g	0,1 g	± 0,1 g	155X155 mm
MINI 5000	5000 g	0,5 g	± 0,5 g	155X155 mm

Optional Accessories:

Cod.	Description
A2-1000	Power supply AC 12V/250mA. (replacement)
A2-1005	4 Rechargeable batteries AA



ESSE3 srl, Via Garibaldi 30
14022 Castelnuovo D.B. (AT)
Tel +39 011 99 27 706
Fax +39 011 99 27 506
e-mail esse3@chierinet.it
web: www.esse3-medical.com

