



**DATA SHEET** 

**Model: ESDS-3019**

**Bluetooth Electronic Stethoscope**



Material	PVC
Length	11.5cm
Color	Black
Size	11.5x1.7x1.4 cm

CE



Bluetooth indicator



Mode indicator

Power/Record




Volume+/Bluetooth

Volume-



Diameter of stainless steel ear hook:  
0.6cm/0.236in


PVC sound guide pipe diameter:  
1.1cm/0.43in



Total length of electronic stethoscope:  
78cm/30.7in

4.5cm/1.77in

Net weight of electronic stethoscope:  
155g (without battery)



App software




Silicone earplug



PVC earplug



Nameplate



Up to 24x amplification  
(at peak frequency, vs, analogic mode)

