



**DATA SHEET** 

**Model: ESF9-0719**

**CPR Mask - pocket resuscitator**



The CPR pocket mask is a non-invasive, latex free mask with inflatable cuff and an integral non-rebreathing valve and efficient filter to prevent direct mouth contact with the patient's face when providing mouth-to-mouth resuscitation. Clear plastic, designed to fit adult or child. Soft inflated cuff. Mask not affected by extreme temperatures. Standard 22 mm. fittings - can be used with bag and mask resuscitator. Flexible dome folds into pocket-size mask. Compact carrying case, with instructions.



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)



CE