

# ESF OPERATING TROLLEY

Constructed from stainless steel 18/10.  
Shelves with up stands on four sides.  
Push handle.  
Mounted on  $\varnothing$  80 mm. casters.  
Knock down construction

Structure en acier inox 18/10.  
Deux plateaux avec 4 bords relevés.  
Poi gnée.  
Monté sur roulettes  $\varnothing$  80 mm.  
Structure démontable.



## MEDICAL TROLLEY / TABLE DE PANSEMENT

CODE	SIZE/DIMENSIONS
<b>ESF2-869/01</b>	<b>70x50x80H cm</b>



## INTRUMENT TROLLEY / GUERIDON PORTE INSTRUMENTS

CODE	SIZE/DIMENSIONS
<b>ESF2-866/01</b>	<b>60x40x80H cm</b>
<b>ESF2-867/01</b>	<b>70x50x80H cm</b>
<b>ESF2-868/0 1</b>	<b>90x60x80H cm</b>



## DRESSING TROLLEY / TABLE DE PANSEMENT

CODE	SIZE/DIMENSIONS
<b>ESF2-1050/01</b>	<b>1 drawer, ss bowl, 70x50x80H cm</b>
<b>ESF2-1051/01</b>	<b>2 drawers, ss bowl, 70x50x80H cm</b>



## DRESSING TROLLEY / TABLE DE PANSEMENT

CODE	SIZE/DIMENSIONS
<b>ESF2-863/01</b>	<b>Rack, wastebin, ss bowl, 90x60x80H cm</b>

# ESF OPERATING ACCESSORIES



## MAYO TABLE / TABLE DE MAYO

CODE	SIZE/DIMENSIONS
<b>ESF2-481/00</b>	<p>Stainless steel construction. Variable height, 4 castors. Swivel casters 50 mm            Dimensions: 45 x 65 x 85/115cm            Capacity: 25Kg</p> <p>En acier inoxydable. Hauteur réglable. Piètement mobile roulettes 50 mm            Dimensions: 45 x 65 x 85/115cm            Capacité: 25Kg</p>



## KICK BUCKET / BAC A DECHETS

CODE	SIZE/DIMENSIONS
<b>ESF2-541/00</b>	<p>Manufactured from 18/10 stainless steel. Mounted on 5 casters Ø 40 mm.</p> <p>-</p> <p>Acier inoxydable 18/10.            Monté sur 5 roulettes Ø 40 mm.</p>



## PLASTER BOWL / CUPULE POUR PLATRE

CODE	SIZE/DIMENSIONS
<b>ESF2-551/00</b>	<p>Manufactured from 18/10 stainless steel w/o welding..</p> <p>-</p> <p>Acier inoxydable 18/10 sans soudure.</p>



## BOWL / CUVETTE

CODE	SIZE/DIMENSIONS
<b>ESF2-535/00</b>	<p>Stainless steel 18/10 Ø 32 cm. Capacity 4 lt.</p> <p>-</p> <p>Acier inoxydable 18/10, Ø 32 cm. Capacité 4 lt.</p>

# ESF OPERATING ACCESSORIES



## BOWL'S STAND

Constructed from chrome plated steel tube.  
Variable height.  
Stainless Steel bowl.

## PORTE CUPULES

Colonne de support en tube d'acier chromé.  
Hauteur réglable.  
Cuvette en acier inoxydable

CODE		
ESF2-522/00	Single basin stainless , steel base	Simple basin , base en acier inoxydable
ESF2-523/00	Single basin , plastic base	Simple basin , base en plastique
ESF2-532/00	double basin stainless , steel base	Double basin , base en acier inoxydable
ESF2-533/00	double basin stainless , plastic base	Double basin , base en plastique



CODE	WASTE BIN Constructed from stainless steel 18/10. Removable inner container of plastic material.	POUBELLE Acier inoxydable 18/10. Seau amovible en plastique.
ESF2-424/00	Capacity 20 lt	Capacité de 20 l
ESF2-425/00	Capacity 14 lt	Capacité de 14 l
ESF2-411/00	Capacity 5 lt	Capacité de 5 l

# ESF OPERATING ACCESSORIES



## **ESF2-490/00 ANAESTHESIA TROLLEY**

Constructed from stainless steel 18/10.  
Upper section with a large open space.  
Six drawers and storage compartment  
completewith shelf and locking door.  
Mounted on Ø 80 mm. casters.

DIMENSIONS : 40x60x110 (H) cm.

## **ESF2-490/00 GUERIDON D'ANESTHESIE**

Structure en acier inoxydable 18/10 comprenant :  
1 étagère, 6 tiroirs, 1 armoire de rangement avec  
serrure.

Monté sur roulettes diam. 80 mm.

DIMENSIONS: 40x60x110 (H) cm.



## **WARD SCREEN**

Chrome-plated steel tubular frame, plastic  
curtains washable and flame retardant.  
Mounted on casters  
Knock down construction.

## **PARAVENT**

Tube d'acier chromé panneaux anti-feu et  
lavables.

Mobile sur roulettes

Construction démontable

CODE	
<b>ESF2-824/05</b>	<b>2 panels , size 100x42x167 cm</b>
<b>ESF2-825/05</b>	<b>3 panels , size 150x42x167 cm</b>
<b>ESF2-826/05</b>	<b>4 panels , size 200x42x167 cm</b>



ESSE3 srl, Via Garibaldi 30  
14022 Castelnuovo D.B. (AT)  
Tel +39 011 99 27 706  
Fax +39 011 99 27 506  
e-mail: [danielec@esse3-medical.com](mailto:danielec@esse3-medical.com)  
web: [www.esse3-medical.com](http://www.esse3-medical.com)

