



# ES-OXY-60SA

Oxygen Concentrator



OPERATING THEATRE




**Application**

- Medical use : connect with ventilator and anaesthesia machine.
- Connecting small central oxygen supply.
- Oxygen supply to clinics and small hospitals with 5-50 beds.
- Optional Booster Pump: 2-8 Bar.

**Technical Specification**

Oxygen Flow: 0-60LPM  
 Oxygen Concentration: 96% ( $\pm 3\%$ )  
 Outlet Pressure: 0.04-0.07Mpa (6-10 PSI)  
 Power: AC110V $\pm 10\%$ ,50/60HZ  
 Sound Level: <60db  
 Power Consumption: < 2200w  
 Net Weight: 160Kgs  
 Size: 860x630x1580

**Schematic**



**ESSE3** srl, Via Garibaldi 30  
 14022 Castelnuovo D.B. (AT)  
 Tel +39 011 99 27 706  
 Fax +39 011 99 27 506  
 e-mail [esse3@chierinet.it](mailto:esse3@chierinet.it)  
 web: [www.esse3.dreamgest.com](http://www.esse3.dreamgest.com)

